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**PLEASE READ THE INSTRUCTIONS CAREFULLY AND THOROUGHLY AND KEEP THEM FOR FUTURE REFERENCE**

## Hip orthosis

### DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG states under its sole responsibility that this product is a Class I medical device that has been manufactured in accordance with the requirements of Directives 93/42/ EEC and 2007/47/EC. These instructions have been written in conformity with these directives, and their aim is to ensure the safe and proper use of this medical device.

### TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

### PRECAUTIONS

It is advised that the pressure exerted by the product does not act on parts of the body suffering from wounds, swelling or lumps. The product should not be tightened too firmly so as to avoid the development of pressure points or compression of the nerves or blood vessels underneath. It is advisable to wear the product over a garment, avoiding direct contact with the skin. If in doubt about how to apply, contact an orthopaedic technician.

Please read the composition of the product on the inside label carefully. The device should not be worn in the vicinity of open flames or strong electromagnetic fields. Wearing a knee-length hose under the Kafo module and a sports shoe with laces is recommended.

### WARNINGS

The product has been designed for the specific indications stated below and must be prescribed by a doctor and fitted by an orthopaedic technician. The orthopaedic technician is responsible for fitting the product and providing information about how to use it safely in accordance with individual requirements.

The product must be fitted with the utmost care to ensure its efficacy, tolerability and correct function. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered.

The manufacturer is no longer liable if the product is used or adapted inappropriately. The orthosis is intended to be used only by a single patient; if it is used for more than one patient, the manufacturer accepts no responsibility or liability in accordance with medical device legislation. Direct contact with the skin may cause reddening or irritation in hypersensitive individuals. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician. The orthopaedic efficacy of the product may only be guaranteed when all of its components are used.

### MAINTENANCE

- Do not bleach  No chemical cleaning
- Do not iron  Do not tumble-dry
- Washing instructions:

- Padding: Wash by hand in lukewarm water (max. 30°C) with neutral soap; rinse thoroughly. Dry away from heat sources
- Rigid parts: Clean with a sponge soaked in lukewarm water (max. 30°C) and neutral soap. Dry with a cloth
- Check at regular intervals to make sure that the screws for flexion/ extension and the abduction setting are tightened
- Do not dispose of the product or any of its components into the environment

### INDICATIONS

- Pre and/or post-surgery (arthroplasty or hip revision)
- Inoperable or insufficiently consolidated femoral neck fractures
- Abduction and flexion-extension control for patients with potential joint problems

### CONTRAINDICATIONS

Currently no known

### CHARACTERISTICS AND MATERIALS

#### Hipo

- Anatomically shaped and perfectly adaptable plastic parts
- Modular system that permits replacement and adaptation of parts
- Available in 4 sizes, right and left Kinesio hinges
- Abduction/adduction adjustable in steps of 7,5°
- Flexion-Extension settings from 0° to 120° in steps of 15°
- The Kinesio Offset articulations are manufactured at a special bending angle to provide more room between the articulation and the hip.
- REF. 2306 only: Using the dual hinge the position of the thigh component can be corrected medially or laterally

### ADAPTING

#### ADAPTING FOR THE DOCTOR / TECHNICIAN

The hip brace is made up of: Hip holder, Kinesio joint brace and thigh holder.

#### 1 Putting on the hip holder:

- Loosen the four screws of the lumbar clips of the hip holder (but do not remove fully).
- Open the front buckle of the hip holder.
- Position the hip holder in such a way that the bracket for the joint brace is placed along the centre line of the patient's hip (fig. A), and then repeat the same process for the hip holder on the opposite side.
- Check to ensure that the hip holder rests on the iliac crest and the brace fits closely to the waist.
- Exert a light pressure on the brace so that it is pushed close to the hip and tighten the four screws of the lumbar clip (this process is best carried out by two people) (fig. B).
- Check to make sure that the connecting strap is exactly centred (this can be carried out using the notches in the connecting strap).
- Set the stomach strap to the appropriate length and hook into place (fig. C). If necessary, cut off any protruding parts and carefully seal the edge using a small flame.

#### 2 Setting the Kinesio joint brace:

- Remove the clear protective cap after undoing the appropriate fixing screw (fig. D).
- Set the flexion-extension as prescribed by the doctor. Undo the screws, set the degree of flexion-extension and tighten the screws again (fig. E).
- Replace the clear protective cap again and secure using the appropriate screw.
- Loosen the screws for the abduction/adduction joint far enough (but do not remove) that both parts of the brace can rotate freely (fig. F).

#### 3 Attaching the joint brace to the hip holder (fig. G):

- Attach using the supplied screws, do not tighten the screws.
- Align the flexion-extension joint in such a way that the centre point is positioned at the head of the femur neck (to align exactly use the Trochanter major as a reference point). Then tighten both fixing screws on the hip holder.

#### 4 Attaching the thigh holder:

- The bracket for the joint brace must be on the side.
- Ensure that the thigh holder is in a medial position on the condyle and the gap proximal to the bar is at least three to four fingers.
- Close the thigh holder using the appropriate Orthèse pour hanche straps. Before doing so, adjust to the correct length (fig. H). If necessary, cut off any protruding parts and carefully seal the edge using a small flame (cigarette lighter).
- Align the bracket for the joint brace in accordance with the joint brace. Secure the brace at the thigh holder using the two supplied screws. Tighten the screws (fig. I).

#### 5 Setting the adduction/abduction (fig. I):

- Bring the extremity to the abduction or adduction position as prescribed by the doctor.
- In order to secure the required position, tighten the screw for the adduction/abduction joint.

**CAUTION:** If the brace is fitted to a bedridden patient, we recommend setting the hip holder as follows:

1. Measure the length of the rear waist circumference (between the centre lines of the hips);
2. Loosen the four screws of the lumbar clips of the hip holder (but do not remove fully). Set the hip holder in such a way that the gap (measured on the exterior circumference) between the bracket for the joint brace and the centre point of the opposite pelvic half is equal to half of the waist circumference of the patient, increased by 3-4 cm. Tighten the screws of the lumbar clips. Use the notches on the connecting strap to centre it.
3. The angling of the lateral cups of the hip holder opposite the connecting strap must be determined by eye until the patient is able to stand up again.

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