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DESIGNED AND DEVELOPED IN SWITZERLAND

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EN – The descriptions and images in this document are for example and commercial purposes only. Orthoservice reserves the right to modify them according to its needs.

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Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland

Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89

info@orthoservice.com · www.orthoservice.com

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Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13

info@orthoservice.de · www.orthoservice.de

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Tel. +48 (0) 34 340 13 10 · (NIP): 94922246785

info@orthoservice.pl · www.orthoservice.pl

OSUK Ltd

Units 47&50 North East BIC

Wearfield

Sunderland Enterprise Park East

Sunderland SR5 2TA

Tel +44 (0) 191 51 66 220

info@orthoservice.uk · www.orthoservice.uk

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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND SAVE THEM

Lower leg extension for hip orthosis Hipo

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Wearing a knee-length hose under the Kafo module and a sports shoe with laces is recommended. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Item	M2320 DX/SX	M2321 DX/SX
Size	Regular	Long
Limb length (groin - heel) cm	63/75	76/93
Colour	white	

Specify right or left

Foot size	S 36/39	M 40/43	L 44/46	XL >46
-----------	---------	---------	---------	--------

MAINTENANCE

- ☒ Do not bleach
- ☒ No chemical cleaning
- ☒ Do not iron
- ☒ Do not tumble-dry

Washing instructions:

- Padding: Wash by hand in lukewarm water (max. 30°C) with neutral soap; rinse thoroughly. Dry away from heat sources
- Rigid parts: Clean with a sponge soaked in lukewarm water (max. 30°C) and neutral soap. Dry with a cloth
- Check at regular intervals to make sure that the screws for flexion/ extension and the abduction setting are tightened
- Do not dispose of the product or any of its components into the environment

INDICATIONS

- Kafo is an additional device which guarantees the control over minimal rotation movements for the hip orthosis Hipo, in this way it ensures that the joint stays in a natural position.

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Improves patient's wearing comfort: lightweight materials, comfortable brace
- The monolateral design reduces the complications of bedsores, especially in geriatric patients
- Maintains the natural position of the limb during flexion of the hip
- Reduces the risk of dislocation
- Guarantees optimum control of flexion-extension, of abduction and internal/external rotation of the hip

ADAPTING

ADAPTING FOR THE DOCTOR / TECHNICIAN

Kafo is a device used to extend the Hipo hip brace. It was designed to create complete stability of the brace and to achieve the appropriate motion control of the lower extremity.

- | | |
|---------------------------------|--|
| 1 - Nylon plate | 9 - Fixed plate |
| 2 - Upper brace (thigh) | 10 - Loop |
| 3 - Knee joint | 11 - Shoe element made from PE |
| 4 - Locking/unlocking mechanism | 12 - Fixing screws at the thigh holder |
| 5 - Lower brace (calf-foot) | 13 - Setting screws |
| 6 - Padding | 14 - Fixing screw |
| 7 - Velcro® straps | 15 - Setting screws |
| 8 - PP holder | 16 - Securing screw |

- Attach the nylon plate to the (previously attached) thigh holder:
 - Undo the two screws that are connecting the joint brace to the thigh holder.
 - Fit the white nylon plate with the supplied fillister head screws 12 to the thigh holder (the screws must secure the nylon plate and the joint brace to the brace). Tighten the screws 13 (fig. A).
- Undo the upper brace from the lower brace (fig. B):
 - Position the two braces at an angle of approximately 90° to one another.
 - Operate the unlocking/locking mechanism in accordance with the arrows, and at the same time, rotate the joint disc (stop) in the UNLOCK direction. Then remove the lower brace from the fork.
- Secure the upper brace to the nylon plate (that was previously secured to the hip holder) using the supplied long screws (13) and the relevant plastic washers. Do not tighten. Align the upper brace along the centre line of the leg (for exact alignment, using the Trochanter major and the Condylus lateralis as reference points). The centre point of the joint must be at the level of the patella axis. Tighten the screws 13 (fig. C).
- Secure the lower brace back to the upper brace. To do this, first loosen the screws 15 (fig. L):
 - Position the two braces at an angle of approximately 90° to one another.
 - Check to make sure that the joint disc is fully rotated in the UNLOCK direction.
 - Insert the lower brace into the fork by operating the locking/unlocking mechanism.
 - Rotate the disc in the LOCK direction until it engages.
- With the extremity stretched, rotate the lower brace until the foot is in the foot section. Then make sure that the heel fits well at the rear edge of the foot section. At the same time, place the plastic shells around the calf and ankle
- Once the brace has been set to the correct length, lock the setting by tightening the screws 15 (fig. N).
- Tighten the straps around the calf and secure with the Velcro® fastener (fig. E). If the straps are too long, cut to length.
- Close the strap at the instep using the Velcro® fastener (fig. F)
- If possible make sure that the alignment has been made correctly so that the joints are not strained incorrectly when walking. If need be, correct the settings.

If the treating doctor has prescribed the locking of the joint at knee height, take the screw from the compartment next to the joint and insert into the securing opening (fig. G). With the brace locked, the flexion is 8°.

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Velcro® è un marchio di fabbrica registrato di Velcro Industries B.V.