OPTION: REF:ICEGEL

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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Patella strap

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/ physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices.In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF.90000	
Size	small	universal
Circumference (cm)	up to 30	from 31
Colour	black	
fit right and left		

MAINTENANCE

⌘ Do not bleach ⌘ No chemical cleaning

⌘ Do not iron ⌘ Do not tumble-dry

⌘ Washing instructions: hand wash at max. 30° C with neutral soap (a sponge is recommended); do not dry in the vicinity of heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Patellar tendinitis and tendinopathy
- Morbus Osgood-Schlatter disease (inflammation of the front tibial apophysis)

CONTRAINDICATIONS

Currently no known

FEATURES AND MATERIALS

- Front elastic strap to adjust pressure and keep it constant also when moving
- Rigid Velcro® surface for stable fixing of the pressure pad or the cooling pad (optional)
- Foam pad with hypoallergenic fabric for easier positioning.
- Skin-friendly gel pad concentrating pressure on tendon, slip resistant
- Velcro® fastening strap, padded for easy application.
- Soft edge for better comfort especially in knee pit
- Optional: pad for cold therapy, with Velcro® fastening system replacing pad

PUTTING ON THE APPLIANCE

- Pull one of the two ends of the strap out of the appropriate loop (fig. A).
- If the pelotte with the pressure pad has not already been positioned, affix this to the securing band using the adhesive points (fig. B).
- Put on the strap in such a way that the pressure pad is located below the knee at the painful position and the cavity of the gel pressure pad is facing upwards (fig. C).
- Feed the free end of the strap (the one you pulled out earlier) back through the appropriate loop (fig. D).
- Also unfasten the end of the Velcro® closure and pull the strap tight evenly (fig. E).
- Secure the ends to the strap using the Velcro® closure (fig. F).
- Check the position of the pressure pad.
- If the strap is uncomfortable, or causes a tingling sensation, loosen the closure a little. If the strap is too loose, it may be tightened up further.

COOLING

- Take the gel pad out of the velour pocket.
- Prior to use, put the gel pad in the freezer for at least two to three hours. (NOTE: DO NOT put the velour pocket in the freezer with the gel pad inserted).
- How to use:
 - Prior to use, put the gel pad into the velour pocket (fig. G). Place the gel pad and the pocket on the Velcro® surface in place of the pressure pad that is used for active usage (fig. H and I)
 - The strap should be applied in exactly the same way as described above
 - Application time: As prescribed by the physician

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Velcro® jest zastrzeżonym znakiem towarowym Velcro Industries B.V.