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 PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Short knee brace for patella instability

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices.In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF.90010					
Size	XS	S	M	L	XL	XXL
Circ. middle of the knee cm	30/33	33/36	36/39	39/43	43/48	48/53
Colour	black					
fit right and left						

MAINTENANCE

- ⌘ Do not bleach
- ⌘ No chemical cleaning
- ⌘ Do not iron
- ⌘ Do not tumble-dry
- ⌘ Washing instructions: hand wash at max. 30° C with neutral soap (a sponge is recommended); do not dry in the vicinity of heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Knee pain due to bruises
- Early stage of arthritis
- Chondropathy/ Chondromalacia patella
- Mild patella instability
- Prevention of patella luxation during sports

CONTRAINDICATIONS

Currently no known.

FEATURES AND MATERIALS

- Multifunctional: reversible for both knees
- C-shaped foam stabilizer for mild instability
- Double-face brace that may be positioned both medial and lateral positions
- Open kneecap
- Black and low profile structure made of orthopedic neoprene
- Comfortable edges

PUTTING ON THE APPLIANCE

DONNING THE SUPPORT FOR THE FIRST TIME WITH A DOCTOR/TECHNICIAN

- Arrange the knee support according to the knee (right/left) and the type of stabilisation (lateral/medial) required. Please note that the RO+TEN logo should always be at the top of the knee support and the C should be on the side that needs to be stabilised; turn the knee support inside out if necessary (fig. C–H).
 - Slide the tubular support onto the knee with the RO+TEN logo facing up (fig. A).
 - Position the knee support so that the kneecap fits perfectly in the opening provided (fig. B).
- ###### SUBSEQUENT WEAR
- Slide the tubular support onto the knee with the RO+TEN logo facing up (fig. A).
 - Position the knee support so that the kneecap fits perfectly in the opening provided (fig. B).

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