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 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

REF. G1104 · G1116 · G1160 · G1163 Aer-OS knee brace

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medcal device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames or strong electromagneticfields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF. G1104 · G1116 · G1160 · G1163					
Size	XS	S	M	L	XL	XXL
Thigh circ. taken 15 cm above mid-patella cm	33/39	39/45	45/52	52/58	58/66	66/75
Length cm REF.G1104 · G1116	34	34	34	34	34	34
Length cm REF.G1160 · G1163	27	27	27	27	27	27
Colour	black					

GenuFIT 04/60/63: fits right and left · GenuFIT 16: specify right or left

MAINTENANCE

 Do not bleach  No chemical cleaning
 Do not iron  Do not tumble dry
 Washing instructions: hand wash at max. 30° C with neutral soap, do not dry in the vicinity of heat sources.
Do not dispose of the product or any of its components into the environment.

INDICATIONS

GenuFIT 04/60/63

- Knee pain due to bruises
- Early stage arthrosis
- Solo GenuFIT 04/63: Chondropathy/ Chondromalacia patella
- Solo GenuFIT 04: Mild patella instability
- Solo GenuFIT 04/63: Prevention of patella luxation during sports
- Solo GenuFIT 60: 1st grade sprains to collateral ligaments
- Solo GenuFIT 60: Tendinitis

GenuFIT 16

- Mid-lateral patella instability
- Patellofemoral arthrosis
- Gonarthitis

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

GenuFIT 04/60/63

- New anatomical, ultra-thin design for a flexible movement and perfect adherence
- Breathable Aer-OS fabric
- New softer Soft-X fabric for kneepit insert (with no central stitch) for better comfort
- Softer and elastic edge
- GenuFIT 04 only: Patella opening with foam stabiliser
- GenuFIT 04 only: Thin lateral stays to keep brace correctly positioned
- GenuFIT 63 only: Reinforced patella stabilizer

GenuFIT 16

- Structure made of Breathable Aer-OS fabric
- Open Kneecap with a lateral J-shaped foam stabilizer
- Medial and lateral spiral stays
- With winches (adjustable)
- Soft-X fabric for kneepit insert for better comfort

PUTTING ON THE APPLIANCE

GenuFIT 60

- 1 Pull the knee bandage over the knee so that the orange label is facing upwards (fig. A).
- 2 Position the knee bandage in such a way that the orange profile is on the back.

GenuFIT 04 · GenuFIT 63

- 1 Pull the knee bandage over the knee so that the orange label is facing upwards (fig. B).
- 2 Position the knee bandage so that the patella cut-out is positioned precisely over the knee cap (fig. C).

GenuFIT 16

- 1 Open the straps on the patella cut-out and provisionally secure the Velcro® closures lightly.
- 2 Pull the knee bandage over the knee so that the orange label is facing upwards (fig. D), and positioned in such a way that the patella cut-out is located precisely above the knee cap and the stabilisation elements are at the sides of the patella cut-out (fig. E). (In the event of medial instability, the stabilisation must be on the medial side of the patella).
- 3 Feed the upper strap through the appropriate buckle, pull tight and secure with the Velcro® closure (fig. F).
- 4 Repeat this process for the lower strap (fig. G).

NOTE: If the hook tape is too long and protrudes over the loop tape, it may be cut to the appropriate length.

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Velcro® è un marchio di fabbrica registrato di Velcro Industries B.V.