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 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

REF.33108-33108A Knee brace with spiral stays
REF.33115-33115A Knee orthosis with hinges

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been designed and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

It is recommended that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagneticfields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never after the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices.In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF.33108-33108A / REF.33115-33115A					
Size	XS	S	M	L	XL	XXL
Thigh circ. taken 15 cm above mid-patella cm	33/39	39/45	45/52	52/58	58/66	66/75
Length cm	34	34	34	34	34	34
Colour	grey					

fits right and left

MAINTENANCE

 Do not bleach  No chemical cleaning
 Do not iron  Do not tumble-dry
 Washing instructions: hand wash at max. 30° C with neutral soap, (REF. 33115 - 33115A only: after removing the hinged stays), do not dry in the vicinity of heat sources.
Do not dispose of the product or any of its components into the environmen.

INDICATIONS

- GenusKILL 08 / 08A only:Early stage knee pain or effects after sprain trauma with mild residual collateral ligament laxity (non-operative knee)
 - GenusKILL 15/15A only: Knee pain or effects after sprain trauma with mild residual collateral ligament laxity (non-operative knee)
 - GenusKILL 08/15 only: Mild mid-lateral patella instability
 - GenusKILL 08/15 only: Patellofemoral arthris
- CONTRAINDICATIONS**
Currently no known

CHARACTERISTICS AND MATERIALS

- Structure made of BreathSkill:
 - weich und breathable
 - 3D weave
 - high resistant mesh
- GenusKILL 08/15 only: FlexSkill front support (bi-layer technical fabric) with special low-profile anti-irritation welding seam
- GenusKILL 08A/15A only: New Velcro® fastening system for easier fastening
- GenusKILL 08/15 only: Patella stabilizer
- CompresSkill support:
 - Stretchy multi-fibre and compression material
 - Ultra-low-profile anti-irritation for better comfort (knee flexion)
 - Fully adaptable
 - Perfect adherence
- GenusKILL 08/08A only: Medial and lateral stays to keep brace correctly positioned.
- GenusKILL 15/15A only: Pair of polycentric hinges for mid-lateral stabilization (flexion limit to 60° or 90°).
- Upper and lower stabilizing straps:
 - Symmetrical fastening system to avoid brace rotation when fastening
 - Articulated central buckle for "automatic" adaptation to leg
 - Rubber tab for easy Velcro® fastening

ADAPTING FOR THE DOCTOR/TECHNICIAN

- Undo all the Velcro® fastenings and fasten them back on themselves for the time being.
- GenusKILL08 - GenusKILL15**
Slide the tubular support onto the knee with the orange label towards the top (fig. A), the kneecap positioned exactly in the middle of the opening, (fig. B) and the stays parallel to each other.
- GenusKILL08A - GenusKILL15A**
Wrap the support around the knee with the orange label facing upwards, taking care to align the opening with the kneecap; fasten the Velcro® straps above and below the knee (fig. C and D).
- GenusKILL15 - GenusKILL15A**
If necessary, remove the knee support, remove the stays from their pockets and adapt them to fit the shape of the leg (fig. E). Next slide the stays back into their pockets and don the knee support once again.
- Pass the two ends of the upper strap through the ring provided, pull symmetrically with both hands, and fasten with the Velcro® (fig. F). Repeat this procedure with the lower strap (fig. G).

PUTTING ON THE APPLIANCE FOR THE PATIENT

- Undo all the Velcro® fastenings and fasten them back on themselves for the time being.
- GenusKILL08 - GenusKILL15**
Slide the tubular support onto the knee with the orange label towards the top (fig. A), the kneecap positioned exactly in the middle of the opening, (fig. B) and the stays parallel to each other.
- GenusKILL08A - GenusKILL15A**
Wrap the support around the knee with the orange label facing up, taking care to align the opening with the kneecap; fasten the Velcro® straps above and below the knee (fig. C and D).
- Pass the two ends of the upper strap through the ring provided, pull symmetrically with both hands, and fasten with the Velcro® (fig. F). Repeat this procedure with the lower strap (fig. G).

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