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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Knee orthosis with hinges and F-E control, 50 cm long wraparound version

DECLARATION OF CONFORMITY
As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS
Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS
We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields.

WARNINGS
It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF.G1150A					
Size	XS	S	M	L	XL	XXL
Thigh circ. taken 15 cm above mid-patella cm	33/39	39/45	45/52	52/58	58/66	66/75
Length cm	50	50	50	50	50	50
Colour	black					

Fits right and left

MAINTENANCE

- ☞ Do not bleach
- ☞ Do not iron
- ☞ Washing instructions: hand wash at max. 30°C with neutral soap. Dry away from heat sources. Do not dispose of the product or any of its components into the environment.
- ☒ No chemical cleaning
- ☒ Do not tumble-dry

INDICATIONS

- Post-surgery for the protection of knee flexion and extension (ligament reconstruction, transposition of the anterior tibial apophysis, fracture of the distal third of femur, fractures of proximal third of tibia, fractures of patella, tendon injuries of extensor system)
- Post-traumatic stage of knee sprains
- Functional treatment of first and second sprain injuries of lateral/medial collateral ligament
- Prevention of recurring sprain LCM and LCL injuries during sports
- Mid-lateral chronic instability
- Instability due to arthrosis. Non-operative knee

CONTRAINDICATIONS
Currently no known

CHARACTERISTICS AND MATERIALS

- Kneecap length: 50cm
- Breathable AirX™ fabric and Velcro® fastening strap
- Safe and easy-to-use snap buckles
- Front opening
- Ultra-low profile for polycentric joint easy to adjust using wedges:
 - Extension: 0°, 10°, 20°, 30°, 45°, 60°
 - Flexion: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°, 105°, 120°
 - Block: 0°
- Opening for rotula
- Softer and elastic edge
- Opening in the kneepit
 - Ultra-low profile anatomic design for a good flexibility and an excellent adherence

PUTTING ON THE APPLIANCE
FIRST APPLICATION

- Open the straps completely:
 - extend them to the maximum length and close them on top
 - open the snap buckles (pressing on the upper and lower grooves) and detach the Velcro® parts positioned near the rod, on the opposite side of the kneecap (fig. A).
- Position the kneecap behind the leg (with the orange label upward), stretch the elastic fabric outward and bring forward the two edges to wrap the leg with extended knee, centring the opening for the rotula as much as possible. Remove the kneecap if necessary, remove the rods from the sheaths, shape and (fig.D). Secure the Velcro® closure tongues (first the upper one [fig.A], then the lower one [fig.B]).
- Remove the kneecap if necessary, remove the rods from the sheaths, shape and reposition them inside the sheaths, and re-apply the kneecap.
- Close all the straps with the snap buckles (fig.D-E, G-H).
- Stretch the straps at the required tension through the specific Velcro®, first those of the thigh (fig.F), then those of the calf (fig. I); close the thigh distal strap first.
- If necessary set the flexion/extension restriction using the wedges (fig. K): Remove the cover, position the wedges at the appropriate degree and close the cover again.
 - Size limits: 0°, 10°, 20°, 30°, 45°, 60°
 - Flex limits: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°, 105°, 120°

SUBSEQUENT APPLICATIONS

- Open the straps completely and detach the Velcro® parts (positioned near the rod) on the opposite side of the kneecap (fig.N).
- Position the kneecap behind the leg (with the orange label upward), stretch the elastic fabric outward and bring forward the two edges to wrap the leg with extended knee, centring the opening for the rotula as much as possible. Secure the Velcro® closure tongues (first the upper one [fig.M], then the lower one [fig.N]).
- Close all the straps with the snap buckles: first those of the thigh (fig.O-P), then those of the calf (fig.Q-R); close the thigh distal strap first.

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