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   **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

**REF. 8104-8105-8106-8107 Knee orthosis for immobilization of 0°**

**REF. 8104F-8105F-8106F-8107F Knee orthosis for immobilization of 20°**

**DECLARATION OF CONFORMITY**

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medcal device.

**TRADEMARKS OF MATERIALS**

Velcro® is a registered trademark of Velcro Industries B.V.

**SAFETY PRECAUTIONS**

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

**WARNINGS**

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient.

To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

**SELECTION/SIZE**

| Item                | REF.8104  | REF.8105  | REF.8106  | REF.8107  |
|---------------------|-----------|-----------|-----------|-----------|
| Immobilisation at   | <b>0°</b> | <b>0°</b> | <b>0°</b> | <b>0°</b> |
| Länge in cm         | 40        | 50        | 60        | 70        |
| Colour              | black     |           |           |           |
| fits right and left |           |           |           |           |

| Item                | REF.8104/F | REF.8105/F | REF.8106/F | REF.8107/F |
|---------------------|------------|------------|------------|------------|
| Immobilisation at   | <b>20°</b> | <b>20°</b> | <b>20°</b> | <b>20°</b> |
| Länge in cm         | 40         | 50         | 60         | 70         |
| Colour              | black      |            |            |            |
| fits right and left |            |            |            |            |

**MAINTENANCE**

⌘ Do not bleach   No chemical cleaning
⌘ Do not iron   Do not tumble-dry
⚠ Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

**INDICATIONS**

- REF. 8104 – 8107:** Post-surgery or post-traumatic immobilization with knee in extension
- REF. 8104/F – 8107/F:** Post-surgery or post-traumatic immobilization with knee in extension of 20°

**CONTRAINDICATIONS**

Currently no known

**CHARACTERISTICS AND MATERIALS**

- Structure made of double-lined soft foam
- Terry cloth interior
- Medial and lateral rigid stays and back stays made of adjustable light alloy (only REF. 8104/F – 8107/F: with flexion of 20°)
- Adjustable medial and lateral panels
- Patella opening
- Washable

**PUTTING ON THE APPLIANCE**

**ADAPTING FOR THE DOCTOR / TECHNICIAN**

- Open the Velcro® tape closures and provisionally secure to each other. Loosen the side supporting plates from the rear plate.
- With the help of the supplied Velcro® tapes, place the rear supporting plate on the centre of the leg (fig. A). One end of the tape should be as close as possible to the edge of the plate.
- Secure the side supporting plates with the Velcro® closure to the rear plate and make sure that the rods are positioned precisely on the sides and the knee cap is located in the cut-out (fig. B). The opening on the front must lie precisely over the knee cap.
- If necessary, bend the rods to correspond to the anatomical profile of the leg.
- Tighten and secure the Velcro® closures using an alternating sequence.

**PUTTING ON THE APPLIANCE FOR THE PATIENT**

- Open the Velcro® tape closures and provisionally secure to each other.
- Close the Velcro® tapes on the rear supporting plate to provisionally secure to the knee immobilisation support. The opening on the front must lie precisely over the knee cap.
- Tighten and secure the Velcro® closures using an alternating sequence (high/caf) starting with the closures closest to the centres of the knee (fig. C).

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