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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Knee orthosis for post-surgery with graduated movement

DECLARATION OF CONFORMITY As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient.

SELECTION/SIZE

Table with columns: Item, Size (Regular, Long, Rehab), Limb length groin - ankle cm, Length rods cm, Colour.

MANTENANCE

- Do not bleach, Do not iron, Washing instructions: Hand wash at max. 30° C with neutral soap. Dry away from heat sources. Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Post-surgery for the protection of knee flexion and extension (ligament reconstruction, transposition of the anterior tibial apophysis, fracture of the distal third of femur, fractures of proximal third of tibia, fractures of patella, tendon injuries of extensor system)

CHARACTERISTICS AND MATERIALS

- Valves with strap-loops and structural hook-and-loop-system, Non-metallic adjusting joint, light and intuitively adjustable, Flexion: from 0° to +120° with intervals of 10°

ADAPTING

- Open all straps and temporarily close them upon themselves.
- Separate the padding from the rods with the straps, closing it at the back with the appropriate Velcro® fasteners.
- Apply the padding such that it fits closely to the leg, closing it at the back with the appropriate Velcro® fasteners.

ADAPTING FOR THE DOCTOR / TECHNICIAN

- Press the orange button and at the same time rotate the transparent cover (fig. 1).
- To adjust flexion (fig. 2): Rotate the opening of the cover to coincide with the flexion stop and remove the pin.

ADAPTING FOR THE DOCTOR / TECHNICIAN

- Remove the rods from the two rods the valves pressing towards the centre, placing the screwdriver in the appropriate notch and pressing outwardly (fig. 8).
- Place the rod on a worktable so that the edge is aligned with the grooves on the rod (which coincides with the desired length). Press downwards till it breaks (fig. 9).

RECOMMENDATIONS FOR THE PATIENT

- Open the upper part (straps and padding) and wash the thigh.
- Close the upper part.
- Open the lower part (straps and padding) and wash the leg.

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