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 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND SAVE THEM**

Graduated movement post-operative knee orthosis, telescopic version

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG states under its sole responsibility that this product is a Class I medical device that has been manufactured in accordance with the requirements of Directives 93/42/EEC and 2007/47/EEC. These instructions have been written in conformity with these directives, and their aim is to ensure the safe and proper use of this medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS

The product has been designed for the specific indications stated below and must be prescribed by a doctor and fitted by an orthopaedic technician.

The orthopaedic technician is responsible for fitting the product and providing information about how to use it safely in accordance with individual requirements. The product must be fitted with the utmost care to ensure its efficacy, tolerability and correct function. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. The manufacturer is no longer liable if the product is used or adapted inappropriately. The orthosis is intended to be used only by a single patient; if it is used for more than one patient, the manufacturer accepts no responsibility or liability in accordance with medical device legislation. Direct contact with the skin may cause reddening or irritation in hypersensitive individuals. In the case of onset of pain, swelling, lumps or other problems, please consult your physician immediately. The orthopaedic efficacy of the product may only be guaranteed when all of its components are used.

SELECTION/SIZE

Item	REF. G1046
Size	one size
Color	black

Fits right and left

MAINTENANCE

⚠ Do not bleach ⚠ No chemical cleaning

⚠ Do not iron ⚠ Do not tumble-dry

🧼 Washing instructions: Rods: wash with a sponge soaked in warm water and mild soap; dry with a cloth.

Padding: remove the padding and wash them in warm water with mild soap; allow to dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Post-operative treatment that requires control of knee flexion-extension ligament reconstruction, anterior tibial apophysis transposition, fractures of the distal third of the femur, fractures of the proximal third of the tibia, fractures of the patella, tendon damage to the extensor apparatus
- Post trauma phase of knee sprains
- Varus or valgus osteotomies

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- The telescopic lateral hinged stay is made from a lightweight aluminium alloy. Telescopic adjustment of the length with easy-to-use buttons on the cuffs.
- Shells with strap-loops and structural auto-fastening system, non-metallic adjusting joint, light and intuitively adjusted:
 - Flexion: from 0° to +120° with intervals of 10°
 - Extension: from -10° to +60° with intervals of 10°
 - Possibility to lock the articulation from 0° to +60° with intervals of 10°
 - Removable and washable condyle protection pads
- Pads in MTP to limit migration of guide and to protect leg from direct contact with the rods; lined only externally, cuttable, for a good fit around the leg

ADAPTING

- Open all straps and temporarily close them upon themselves.
- Separate the padding from the shell with the straps.
- Follow these steps to adjust the length of the telescopic stays for thigh and calf:
 - Raise the button provided on each cuff (fig. A1)
 - Slide the cuff until it is the desired length while keeping the button raised (fig. A2)
 - Release the button so that it returns inside the corresponding opening on the metal part.
- Apply the padding such that it fits closely to the leg, closing it at the back with the appropriate Velcro® fasteners: first the thigh padding, then the calf padding (fig. B). In order to position the padding at the optimal height, place the orthosis next to the leg. If necessary the padding may be shortened to adapt it to the patient's legs (shorten the padding along the thinner parts in such a way that the padding closes properly).
- Model the rods (by bending them) to fit them to the shape of the patients' legs (fig. C).
- Adjust the flexion extension as described (points 1-5, hinge adjustment):perform the same operation on both joints.
- Fix the rods to the pads (fig. D): the rods should follow the median line of the leg and the joints should be aligned with the patella so as to rest on the condyles. To align the rods with the median line of the leg, use the lateral malleolus and the large trochanter as reference points.
- Ensure that the Velcro® straps are perfectly aligned with the corresponding padding.
- Ensure that the two joints are at the same height.
- Velcro® fasten the straps after inserting them in the appropriate rings (fig. E-F): fasten first the thigh straps beginning with the one closest to the knee; then proceed by fastening the calf straps. Tighten the straps so as to avoid vertical migration of the guide.

NOTE: The bending procedure should be performed by hand, with the help of the grips or laying the orthosis on a rounded body (chair, edge of a bed); do not use your own knee as a lever.

ADAPTING FOR THE DOCTOR / TECHNICIAN HINGE ADJUSTMENT

- Press the orange button and at the same time rotate the transparent cover (fig. 1).
- Pour ajuster la flexion (fig. 2):
 - Rotate the opening of the cover to coincide with the flexion stop and remove the pin
 - Rotate the cover and the pivot to the hole with the degree of flexion required (shown on the mechanism) and insert the pin
- Rotate the cover again (fig. 3).
- To adjust extension (fig. 4):
 - Rotate the opening of the cover to coincide with the flexion stop and remove the pin
 - Rotate the cover and the pivot to the hole with the degree of extension required (shown on the mechanism) and insert the pin
- Rotate the cover until the orange safety button clicks (fig. 5).

RECOMMENDATIONS FOR THE PATIENT

To facilitate washing the immobilized limb without removing the brace, proceed as follows:
1 Open the upper part (straps and padding) and wash the thigh
2 Close the upper part
3 Open the lower part (straps and padding) and wash the leg
4 Close the lower part
5 Tighten the strap to prevent vertical migration of the brace.
To wash the limb we recommend using a damp sponge.

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