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AIR PAD

SINGLE PATIENT
MULTIPLE USE

 PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Knee orthosis for post-surgery with graduated movement, PCL version

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagneticfields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices.In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF.G1044
Size	one size
Stay length cm	56
Colour	black

Fits right and left

MAINTENANCE

- | | |
|---|--|
|  Do not bleach | <input checked="" type="checkbox"/> No chemical cleaning |
|  Do not iron | <input checked="" type="checkbox"/> Do not tumble-dry |
|  Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources. Stays: clean with a sponge dipped in lukewarm water and mild soap; dry with a cloth. Pads: remove the pads and wash them in lukewarm water with mild soap; leave to dry away from direct heat. | |
- Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Post-surgery when knee flexion and extension need to be limited and controlled (PLC ligament reconstruction)
- Post-traumatic stage of knee sprains

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Inflatable retro-tibial pad to relieve pressure on PCL
- Non-metallic adjusting joint, light and intuitively adjusted
- Flexion: from 0° to +120° with intervals of 10°
- Extension: from -10° to +60° with intervals of 10°
- Possibility to lock joint from 0° to +60° with intervals of 10° structural hook-and-loop system (not bonded Velcro®)
- Frame of hinges made of light aluminum alloy
- MTP Pads to limit brace movement and protect leg from direct contact with hinge; lined only externally, cuttable and perfectly adaptable to leg circumference
- Front and back straps to ensure the stays position along the median line of leg
- Removable and washable condyle protection pads

ADAPTING FOR THE DOCTOR / TECHNICIAN

- Undo all the straps and fasten them back on themselves for the moment.
- Separate the padding from the stays with the straps.
- Check that the retro-tibial pad is completely deflated.
- Attach the pads securely to the leg, fastening them at the front with the Velcro straps provided: first fasten the thigh straps and then the calf straps (fig. A). Place the stays next to the leg to make it easier to position the bands at the right height. If necessary, trim the pads to adapt to the circumference of the leg (cut along the lines where the pads meet so that they do not open up).
- Bend the stays so they can be adapted to fit the shape of the patient's leg (fig. B).
- Adjust the flexion and extension, as described (points 1–5, adjusting the hinges): use the same procedure for both hinges.
- Attach the stays to the pads (fig. C): the stays should follow the median line of the leg and the hinges should be aligned with the kneecap so that they rest on the condyles. Use the lateral malleolus and the greater trochanter as reference points to help align the stays with the median line of the leg.
- Check that the straps are perfectly matched up with the pads.
- Check that the two hinges are at the same height.
- Pass the Velcro® straps through the rings provided and fasten: first fasten the thigh (fig. D) straps beginning with the one nearest the knee, and then fasten the calf straps (fig. E). The straps should be tightened enough to prevent the brace moving up or down.
- Inflate the retro-tibial pad as directed using the pump provided (fig. F). Warning: over inflating the pad can block the circulation.
NOTE: the bending procedure should be performed by hand, with the help of the grips or laying the orthosis on a rounded body (chair, edge of a bed) ; do not use your own knee as a lever.

HINGE ADJUSTMENT

- Press the orange button and at the same time rotate the transparent cover (fig. G).
- To adjust flexion (fig. H):
 - Rotate the opening of the cover to coincide with the flexion stop and remove the pin
 - Rotate the cover and the pivot to the hole with the degree of flexion required (shown on the mechanism) and insert the pin
- Rotate the cover again (fig. I).
- To adjust extension (fig. L):
 - Rotate the opening of the cover to coincide with the extension stop and remove the pin
 - Rotate the cover and the pivot to the hole with the degree of extension required (shown on the mechanism) and insert the pin.
- Rotate the cover until the orange safety button clicks (fig. M).

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