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EN, FR PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Orthosis for gonarthrosis

DECLARATION OF CONFORMITY As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or itches. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Table with 2 columns: Attribute (Item, Size, Colour) and Value (REF.91014 right or left, one size, grey)

RIGHT = for medial gonarthrosis of the right knee (varus right knee) or for lateral gonarthrosis of the left knee (valgus left knee). LEFT = for medial gonarthrosis of the left knee (varus left knee) or for lateral gonarthrosis of the right knee (valgus right knee).

MAINTENANCE

- ☒ Do not bleach ☒ No chemical cleaning ☒ Do not iron ☒ Do not tumble-dry ☒ Washing instructions:

Rams: Clean with a sponge dipped in lukewarm water and mild soap. Dry with a cloth. Pads: Remove the pads and wash them in lukewarm water with mild soap. Leave to dry away from direct heat. Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Moderate unicompartamental gonarthrosis. • Post-surgery cartilage or meniscus treatments.

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Light structure made of foam and aluminum alloy • Single-sided hinge for improved and effective brace application. Lever mechanism to relieve mid-lateral knee • Hoop-and-lock straps to apply the brace without using Velcro® straps (first setting) • Soft, removable and washable pads • Included accessories: - Additional condyle pad; - wedges to limit extension of 0°, 5°, 10°, 15°, 20° or 25° and flexion to 45°, 60°, 75° or 90°

ADAPTING

INITIAL FITTING

- Open all straps: release the locking levers using the appropriate button (1a+1b). Check that the Velcro® straps are set to maximum length.
2 Sit the patient with the knee flexed at approximately 45° and position the brace on the opposite side of the knee that requires the low relieved; check that the hooks on the levers face forwards and move the centre of the joint so that it is about 1-2 cm above the middle of the kneecap (fig. 2). The upper strap should be number 4.
3 Fasten the straps:
- Attach the levers to the corresponding hooks on the brace and fasten them by pushing until they click the close button (fig. 3); start with the number 1 strap (cafl, proximal) and then fasten the others following the numbers on the buckles.
- Adjust the length of the Velcro® straps using the rubber tips (fig. 4) starting with the number 1 strap (cafl, proximal) and then adjust the number sequence.
- If the strap is still too long, it can be adjusted further inside the brace. Remove the strap padding, undo the strap and adjust the length; if needed, the hook Velcro® can be removed and the excess cut. Attach the end Velcro® part and reattach the padding.
4 Insert the key supplied in the appropriate mechanisms and rotate it until the desired medial / lateral thrust is obtained (fig. 5a + 5b); it should be adjusted by a ¼ turn each time.
5 Ask the patient to walk. Check that the condyle pad is at the same height as the kneecap, and that the knee brace is correctly positioned and feels comfortable. Adjust the brace further if necessary.

OPTIONAL ACCESSORIES

- An optional additional condyle pad is included in the packaging.
7 The brace also comes with wedges to limit (if prescribed) the degree of extension to 5°, 10°, 15°, 20° or 25° and flexion to 45°, 60°, 75° or 90°.
- Remove the screws that hold the flexion/extension wedges in place and remove the wedges if they have already been fitted (the 0° extension wedge is pre-mounted on the brace).
- Insert the wedges with the desired degree of flexion/extension so that the notches on the wedges match up with the corresponding points on the knee brace. Replace the screws and tighten them fully.

SUBSEQUENT WEAR

- Open all straps: release the locking levers using the appropriate button (fig.1a+1b).
2 Sit down with the knee flexed at approximately 45°, position the brace on the opposite side of the knee to where pain is felt with the lever hooks facing forwards, move the centre of the joint so that it is about 1-2 cm above the middle of the kneecap (fig. 2). The upper strap should be number 4.
3 Attach the levers to the corresponding hooks on the brace and fasten them by pushing until they click of the close button (fig. 3); start with the number 1 straps (cafl, alto) and then fasten the others following the numerical sequence.

CAUTION: do not undo the Velcro® straps (which were adjusted during the initial fitting), instead use only the levers to fasten the brace.

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