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 PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Neoprene ankle support with spiral stays

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.
Breathe-o-prene® is a registered trademark of Nexus Foams Ltd.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use. Individuals who are suffering from dermatitis or allergies to products containing neoprene should not use this product.

SELECTION/SIZE

Code	REF.8280			
Size	S	M	L	XL
Ankle circumf. cm	15/20	20/25	25/30	30/36
Colour	blue			

Fits right and left

MAINTENANCE

- ⌘ Do not bleach
 - ⌘ No chemical cleaning
 - ⌘ Do not iron
 - ⌘ Do not tumble-dry
- 🧼 Washing instructions: hand wash at max. 30° C with neutral soap (a sponge is recommended); do not dry in the vicinity of heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Minor tibiotarsal sprains 1st degree
- Moderate inflammatory conditions (peroneal tendinopathy)
- Mild joint pain
- Tibiotarsal arthrosis

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Double-lined orthopedic neoprene
- Spiral stays for support on both sides
- Fastens around the ankle with Velcro® straps
- Open heel
- Washable

PUTTING ON THE APPLIANCE

- 1 Open all straps and temporarily close them back on themselves.
- 2 Put the anklet with the inscription on the foot and the velcro straps on the leg.
- 3 Position the anklet so that the heel fits snugly into the hole (A).
- 4 Tighten and close the Velcro straps, the lower one first (B) and then the upper one (C)

REF.8280

Neoprene ankle support with spiral stays

object 80

ORTHOSERVICE
RO+TEN
Take care feel better

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**Headquarter: ORTHOSERVICE AG**
**ORTHOSERVICE DEUTSCHLAND GmbH**
RO+TEN S.r.l.