malleoLITE

REF. 33080

Knöchelbandage mit elastischen Gurten
Ankle support with elasticated figure-8 straps
Chevillière avec sangles élastiques en 8
Бандаж голеностопного сустава с эластичными затяжками „восьмеркой“
Cavigliera con tiranti elastici a 8

ORTHO SERVICE
R+TEN

Headquarter: ORTHOSERVICE AG

Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland

Tel. 0041 91 8220088 · Fax 0041 91 8220089

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ENGLISH

Please read the instructions carefully and thoroughly

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares with sole responsibility that this product is a Class I medical product and has been manufactured in accordance with the basic requirements of the directives 93/42/EEC and 2007/47/EEC. These operating instructions were created in accordance with the specification of the above-mentioned directives. They aim to ensure the safe and correct use of the medical product.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SELECTION/SIZE

REF. 33080	S	M	L
Size			
Circ. of ankle from base of heel across the malleoli cm	24-28	28-33	33-38
Colour:	black / Fits right and left		

MAINTENANCE

Do not bleach No chemical cleaning

Do not iron Do not tumble-dry

Washing instructions: Hand wash in lukewarm water (max 30°C) with mild soap; rinse thoroughly.
Leave to dry away from direct heat.

Dispose of this product and its parts responsibly.

PRECAUTIONS

It is advised that the pressure exerted by the product does not act on parts of the body suffering from wounds, swelling or lumps. The product should not be tightened too firmly so as to avoid the development of pressure points or compression of the nerves or blood vessels underneath. It is advisable to wear the product over a garment, avoiding direct contact with the skin. If in doubt about how to apply, contact an orthopaedic technician. The device should not be worn in the vicinity of open flames.

WARNINGS

The product may only be used in the event of medical indication or when prescribed by a physician and must be applied by an orthopaedic technician, the competent authority figure for the application as well as for information on the safe use, according to individual needs. To ensure the effectiveness, tolerability, and correct operation, it is necessary that the device is put on with the utmost care. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. Use of the product is recommended for a single patient. With hypersensitive patients, direct contact with skin may cause redness or irritation. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician.

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ENGLISH

FEATURES AND MATERIALS

- Made with Sensitive® fabric which is stretchy, breathable, ladder-proof and antistatic.
- The anatomic, open design is ever-so easy to put on.
- Fastens at the back with Velcro®.
- Elasticated crossover instep straps.

INDICATIONS

- Grade 1 ankle sprains.
- Tendinopathy of the instep and peroneal tendons.
- Mild joint pain.

CONTRAINDICATIONS

There are no known contraindications at this time.

DONNING THE SUPPORT

- 1 Undo all the Velcro® fastenings and slide the forefoot into the ankle support (fig. A).
- 2 Fasten the Velcro® at the back so that the ankle support fits snugly (fig. B).
- 3 Grasp the lateral elasticated strap, pull it over the top of the foot, and fasten it on the medial side of the foot with the Velcro® (fig. C).
- 4 Grasp the medial elasticated strap, pull it over the top of the foot, and fasten it at the side with the Velcro® (fig. D).

WARRANTY

This warranty is governed by the statutory provisions of the country in which the product was purchased. Where relevant that country's regulations concerning warranties and the relationship between the seller and buyer shall apply. If you believe that the terms of the warranty may be applicable, first contact the retailer where the product was purchased. Orthoservice will repair or replace all or part of the device and its accessories if defects in materials and workmanship arise during the period of six months from the date of purchase.

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