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 PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Ankle support with elasticated figure of 8 straps

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF.33080
Size	S M L
Circ. of ankle from the malleoli cm	24/28 28/33 33/38
Colour	black
fits right and left	

MAINTENANCE

- ☒ Do not bleach ☒ No chemical cleaning
- ☒ Do not iron ☒ Do not tumble-dry
- ☒ Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- 1st degree tibiotarsal sprains
- Peroneal or instep tendinopathy
- Mild joint pain.

CONTRAINDICATIONS

There are no known contraindications at this time.

FEATURES AND MATERIALS

- Made of Sensitive® fabric which is stretchy, breathable, rip-proof and antistatic.
- The anatomic and open design make it extremely easy to put on.
- Fastens at the back with Velcro®.
- Elasticated straps cross over the instep.

DONNING THE SUPPORT

- 1 Undo all the Velcro® fastenings and slide the forefoot into the ankle support (fig. A).
- 2 Fasten the Velcro® at the back so that the ankle support fits snugly (fig. B).
- 3 Grasp the lateral elasticated strap, pull it over the top of the foot, and fasten it on the medial side of the foot with the Velcro® (fig. C).
- 4 Grasp the medial elasticated strap, pull it over the top of the foot, and fasten it at the side with the Velcro® (fig. D).

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