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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

### Neoprene ankle support with elasticated figure of 8 and spiral stays

#### DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

#### TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.  
Breath-o-prene® is a registered trademark of Nexus Foams Ltd.

#### SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

#### WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

#### SELECTION/SIZE

Code	REF.8281			
Size	S	M	L	XL
Circ. instep to heel cm	26/28	29/31	32/34	35/37
Colour	black			

Fits right and left

#### MAINTENANCE

- Do not bleach
- No chemical cleaning
- Do not iron
- Do not tumble-dry
- Washing instructions: hand wash at max. 30° C with neutral soap (a sponge is recommended); do not dry in the vicinity of heat sources.

Do not dispose of the product or any of its components into the environment.

#### INDICATIONS

- Conservative treatment following grade 2 sprain trauma
- Acute tendinopathies
- Postoperative rehabilitation
- Aid to resuming sporting activity after grade 2 and 3 sprain trauma
- Chronic joint pain caused by tibiotarsal arthrosis

#### CONTRAINDICATIONS

Currently no known

#### CHARACTERISTICS AND MATERIALS

- Double-lined orthopaedic neoprene
- Spiral reinforcing rods on both sides
- Open heel
- Crossover instep straps with Velcro® fastenings pass under the sole of the foot and around the malleoli for maximum stability
- Fastens at the top with a Velcro® strap

#### PUTTING ON THE APPLIANCE

- 1 Open all Velcro® closures and temporarily close them slightly.
- 2 Pull the bandage over the foot so that the heel is exactly in the middle of the heel opening, and subsequently close the Velcro® closure strap (fig. A).
- 3 Positioning of the lateral Velcro® closure:
  - grip the sewn Velcro® strap laterally, guide over the instep and pull downwards to the opposite side (fig. B);
  - guide under the sole of the foot before the heel and close vertically above the outer ankle (fig. C).
- 4 Positioning of the Velcro® closure on the inside:
  - grip the sewn Velcro® strap medially, guide over the instep and pull downwards to the opposite side (fig. D);
  - guide under the sole of the foot before the heel, tension and close vertically above the inner ankle (fig. E).
- 5 Pass the upper strap through the buckle provided, fasten with the Velcro® (fig. F).

**NOTE:** The Velcro® straps can be tensioned to different extents, depending on which degree of stabilisation is required.

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