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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Bivalve ankle brace with foam padding and latero-medial strap

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF.A1028				
Size	XS	S	M	L	XL
Shoe size	30/34	35/38	39/42	43/47	> 47
Colour	black				

specify right or left

MAINTENANCE

- Do not bleach
- No chemical cleaning
- Do not iron
- Do not tumble dry
- Washing instructions: hand wash at max. 30°C with neutral soap (a sponge is recommended); do not dry in the vicinity of heat sources

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Grade 2 and 3 sprains
- Post-surgical functional recovery
- Functional recovery following malleolar fractures
- Chronic instability of the tibiotarsal and in sporting activity

KONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Thin Breath-o-prene® structure
- Medial and lateral plastic valves are inserted in the brace and are protected by soft MTP padding
- Additional rear padding for improved comfort
- Rear elastic fastening strap
- Latero-medial reinforcing and fastening strap (can be shortened to fit)

PUTTING ON THE APPLIANCE

- 1 Open all the Velcro® straps and insert the forefoot in the ankle brace (fig. A).
- 2 Align the shells on the sides of the ankle and fasten the lower rear strap (fig. B).
- 3 Tighten the diagonal strap at the front (fig. C), pass it around the brace (fig. D), then fasten it with the Velcro® strap so that the ankle brace fits well (fig. E).

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Velcro® jest zastrzeżonym znakiem towarowym Velcro Industries B.V.
Breathe-o-prene® jest zastrzeżonym znakiem towarowym Nexus Foams Ltd.