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**PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

**Bivalve ankle brace with inflatable air padding**

**DECLARATION OF CONFORMITY**

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

**TRADEMARKS OF MATERIALS**

Velcro® is a registered trademark of Velcro Industries B.V.

**SAFETY PRECAUTIONS**

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin.

If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician.

Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

**WARNINGS**

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient.

To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation.

The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices.In hypersensitive individuals, direct contact with the skin may cause redness and irritation.

In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

**SELECTION/SIZE**

<b>Code</b>	<b>REF. A1029</b>
<b>Size</b>	<b>Universal</b>
<b>Colour</b>	white

Fits right and left

**MAINTENANCE**

⚠ Do not bleach  No chemical cleaning

⚠ Do not iron  Do not tumble-dry

⚠ Washing instructions: hand wash at max. 30° C with neutral soap (a sponge is recommended); do not dry in the vicinity of heat sources (or, better, dry with a cloth).

Do not dispose of the product or any of its components into the environment.

**INDICATIONS**

- Grade 2 and 3 sprains
- Post-surgical functional recovery
- Functional recovery following malleolar fractures
- Chronic instability of the tibiotarsal and in sporting activity

**CONTRAINDICATIONS**

Currently no known

**CHARACTERISTICS AND MATERIALS**

- Inflatable air pads
- Anatomically shaped lightweight durable plastic bivalve structure
- Velcro® adjustment to fit heel width
- Easy to wear thanks to full opening
- Pump included

**PUTTING ON THE APPLIANCE**

- Open all Velcro® straps.
- Ensure that the inner inflatable air pads are positioned correctly on the plastic structure.
- Position the heel element under the heel in such a way that the rounded section of the heel element (fig. A) is pointing to the front of the foot.
- Position the support shells at the sides of the ankle bone; if necessary adjust the width of the heel element using the Velcro® closure designed for this purpose, and adjust the width of the heel (fig. B).
- Close the orthosis by pulling firmly on the Velcro® straps:
  - First, secure the Velcro® strap near the ankle bone (fig. C)
  - Then close the upper Velcro® strap (fig. D).
- If necessary inflate or deflate air pad (using the supplied pump designed for this purpose) until the maximum wearing comfort and best hold has been achieved.

**Note:** Do not pump the pad up too much as this may cause circulation problems.

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