malleonite

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 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

Night positioning ankle brace

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. The device should not be worn in the vicinity of open flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF.90062	
Size	S/M	L/XL
Shoe size	36/42	43/48
Color	black	

Right and left

MAINTENANCE

-  Do not bleach
-  Do not iron
-  No chemical cleaning
-  Do not tumble-dry
-  Washing instructions: hand wash at max. 30° C with neutral soap.

Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Plantar fasciitis
- Achilles tendinopathy
- Foot drop
- Pain upon awakening

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- The hinged polymer splint can be positioned for personalised support using the adjustable plantar flexion strap
- Comfortable, soft padding
- The adjustable Velcro® calf and foot straps are padded for optimum comfort.

PUTTING ON THE APPLIANCE

- Undo all the straps and fasten them back on themselves for the moment.
- Slide the foot into position and rest the splint of the top of the foot and the ankle (fig. A).
- Pull the Velcro® calf strap and fasten (fig. B).
- Next pull the foot (fig. C) strap and fasten.
- [For the first time of wearing only] Adjust the dorsal flexion using the Velcro® strap (fig. D)
- If greater stabilisation is required, an additional strap is provided that can be placed around the ankle (fig. E).

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