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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Walker REF.400 - REF.450

DECLARATION OF CONFORMITY

ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. It is advisable to wear the product over a garment, avoiding direct contact with the skin. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF.400 - REF.450
Size	XS S M L XL
Shoe size	33/36 36/39 39/44 44/47 > 47
Colour	black

Fits right and left

Option: REF.501 Wedges for left and right heels.

Option: REF. 501

MAINTENANCE

- ☒ Do not bleach ☒ No chemical cleaning
- ☒ Do not iron ☒ Do not tumble-dry

Washing instructions:

- Padding: hand wash with mild soap in warm water (max. 30°C), and rinse thoroughly. Leave to dry away from sources of heat.
- Rigid components: wipe with a sponge dipped in warm water (max. 30°C) with mild soap. Dry with a cloth.

Do not dispose of the product or any of its components into the environment.

INDICATIONS REF.400

- Malleolar fractures conservative treatment
- Grade 3 sprains
- Post-surgical care for osteosynthesis of malleolar fractures, ankle ligament recontruction and Achilles tendon tenorraphy
- Stable fractures of the foot, ankle, tibia and fibula distal

INDICATIONS REF.450

- Soft tissue injuries
- Mild to severe ankle sprains
- Injuries and pain in the forefoot, foot, and heel
- Stable fractures in the foot and/or ankle
- Post-operative care

CONTRAINDICATIONS

The application of this brace is advised only when it is documented that the fracture is stable and there are acceptable limits of angular and rotational misalignments.

CHARACTERISTICS AND MATERIALS

- Structure made from lightweight, durable plastic
- Non-slip rubber sole
- Inner padding of the sole is made from non-toxic Evazote
- Struts held in place with Velcro® for a firm grip on the padding
- Soft foam pad with double lining in brushed nylon for Velcro® grip
- Velcro® fastening straps
- Washable padding
- Available in high (REF.400) and low (REF. 450) version

PUTTING ON THE APPLIANCE

- Undo all the straps and fasten them back on themselves for the moment. Remove the padding from the foot. Slide your foot into the padding so that the back fits snugly around your heel.
- Fasten the padding and Velcro® over the foot and then around the leg (fig. A). The entire length of the padding should fit snugly, but not so tight that it restricts circulation.
- Open the struts of the boot with both hands and slip the boot on over the padding, aligning the struts with the central line of the ankle (fig. B). Extra padding can be added if necessary to ensure a better and more comfortable fit. The opening on the front must lie precisely over the knee cap.
- Remove the separators marked with "PULL UP" and pull them upwards so that the padding can be fastened with Velcro® to the rigid struts (fig. C).
- Fasten the boot with the straps; start with the strap over your foot and work your way up the leg (fig. D-E).

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