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SINGLE PATIENT
MULTIPLE USE

RADIO
LUCENT

 PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Rigid walker boot (open heel)

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF. MM100		
	S	M	L
Size			
Shoe size	< 38	39/42	43/47
Colour	black		

Fits right and left

MAINTENANCE

-  Do not bleach No chemical cleaning
-  Do not iron Do not tumble-dry

 Washing instructions: hand wash at max. 30°C with neutral soap.

Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Conservative treatment of fractures to the ankle
- 3rd degree distortion traumas of the ankle
- Post-operative rehabilitation of ankle fractures, stabilisation of ankle joint ligament injuries, Achilles tendon injuries
- Stable foot, ankle, tibia and fibular distal fractures

KONTRAINDICATIONS

The application of this protector is advised only when it is documented that the fracture is stable and there are acceptable limits of angular and rotational misalignments.

CHARACTERISTICS AND MATERIALS

- Structure made from lightweight, durable plastic; one piece design
- Non-slip rubber sole
- Struts held in place with Velcro® for a firm grip on the padding
- Open on the heel
- Soft foam pad with double lining in brushed nylon for Velcro® grip
- Velcro® fastening straps
- Washable padding

PUTTING ON THE APPLIANCE

- 1 Loosen the straps and remove the padding inside the boot
- 2 Insert your foot into the padding, so that the back of the pad is adherent to the heel. Attach and wrap the foot straps, then the part of the padding relating to the leg. The padding must be snug throughout its entire length length, but must not impede circulation.
- 3 Open the uprights with both hands and slip the boot, aligning the uprights with the central axis of the ankle.
- 4 Attach the boot straps, starting from the toe and working up the leg.

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