 LÖDZU, UZMANIGI IZLASIET UN SAGLABÄJİET INSTRUKCIJAS

REF. 500/510 (short) Pastaigu zābaks ar piepūšamu gaisa polsterējumu (airstep walker - lat.)

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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

REF. 500/510 (short) Walker with inflatable air padding

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Item	REF.500/510(short)		
Size	S	M	L
Shoe size	35/39	39/44	>44
Colour	grey		
fits right and left			

MAINTENANCE

- ☒ Do not bleach
- ☒ No chemical cleaning
- ☒ Do not iron
- ☒ Do not tumble-dry
- ☒ One patient only
- ☒ Washing instructions:
- Padding: Wash by hand in lukewarm water (max. 30°C) with neutral soap; rinse thoroughly. Dry away from heat sources
- Rigid parts: Clean with a sponge soaked in lukewarm water (max. 30°C) and neutral soap. Dry with a cloth

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Conservative treatment of malleolar fractures
- Severe sprains (grade 3)
- Post-surgical phase of osteosynthesis of malleolar fractures
- Ligamentous reconstructions of the ankle and of the tenorhaphy of the Achilles tendon
- Stable fractures of the foot, ankle, tibia and fibula distal

CONTRAINDICATIONS

The application of this brace is advised only when it is documented that the fracture is stable and there are acceptable limits of angular and rotational misalignments

CHARACTERISTICS AND MATERIALS

- Reinforcing shells on the front and back, for increased stability
- Inflatable air pads for optimum fit to the limb and improved stability
- Structure in lightweight, sturdy plastic

ADAPTING FOR THE DOCTOR / TECHNICIAN

- 1 Open all the straps and spread them apart (on the back). Fasten them temporarily on themselves. Remove the padding inside the plastic boot and open it.
- 2 Put on the padding so that the heel rests firmly up against the rear section and then close using the Velcro® fastening (Fig. A).
- 3 Put the plastic boot on the leg from the bottom, spreading the plastic splint (Fig. B).
- 4 Position it so that the splint is vertical and centered on the ankle (Fig. C).
- 5 Close the Velcro® straps, starting from those closest to the toes (Fig.D) and moving upward (Fig. E).
- 6 Inflate the air pads (Fig. F and G) and ensure comfortable fit. To inflate: Insert the small pump (Fig. H) with the opening IN into the desired valve and inject air by using the pump. To deflate: Insert the small pump with the opening OUT into the desired valve. It is then possible to use the pump actively to pump the air out or leave it in place and press on the padding until it is completely deflated. The healthcare personnel assigned should teach the patient how to obtain adherence, comfort and adequate compression, for when the boot is reapplied later by the patient.
- 7 If necessary, use the additional padding supplied to adapt the splint to the patient's leg (Fig. I).

Important: Remove the boot only when prescribed by the healthcare operator (to remove the boot follow the above sequence in reverse).

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