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RADIO
LUCENT

AIR PAD

airstep tight walker
diabetic
airstep tight walker
diabetic

REF. 530

Walker with inflatable air padding and Diabetic Set

REF. 520

Walker with two inflatable air padding and Diabetic Set

REF. 505

Diabetic set

for Walker

**ORTHOSERVICE
RO+TEN**
Take care feel better

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OSUK Ltd
Units 47&50 North East BIC
Wearfield
Sunderland Enterprise Park East
Sunderland SR5 2TA
Tel +44 (0) 191 51 66 220
info@orthoservice.uk · www.orthoservice.uk

RO+TEN S.r.l.
Sede legale: Via Marco De Marchi 7 · I-20121 Milano (MI) · Italia
Sede operativa e amministrativa:
Via Comasina 111 · I-20843 Verano Brianza (MB) · Italia
Tel. 0039 0362 18 51 200
info@roplusten.com · www.roplusten.com

**ORTHOSERVICE
RO+TEN**
Take care feel better

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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

REF. 530 - REF. 520: Walker with (two) inflatable air padding
REF. 505: Diabetic Set

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Item	REF.520 - REF.530 - REF.505		
Size	S	M	L
Shoe size	35/39	39/44	>44
Colour	grey		

fits right and left

MAINTENANCE

- Do not bleach
- Do not iron
- One patient only
- No chemical cleaning
- Do not tumble-dry
- Washing instructions:
- Padding: Wash by hand in lukewarm water (max. 30°C) with neutral soap; rinse thoroughly. Dry away from heat sources
- Rigid parts: Clean with a sponge soaked in lukewarm water (max. 30°C) and neutral soap. Dry with a cloth

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Particularly suited for use with patients with sensitivities, for example, diabetes mellitus and the associated diabetic foot syndrome
- Ulcer on the metatarsus, malleolus, hallux or heel
- Charcot foot

CONTRAINDICATIONS

REF.505: No known

REF.520-530: The application of this brace is advised only when it is documented that the fracture is stable and there are acceptable limits of angular and rotational misalignments

CHARACTERISTICS AND MATERIALS

REF.505

- Adaptable foam inserts for partial pressure relief of the foot sole
- Special stockings for diabetes mellitus and the associated diabetic foot syndrome
- The diabetic set is designed for use with the Airstep Walker and Airstep tight Walker products

REF.520

- Adaptable foam inserts for partial pressure relief from the sole
- Special stockings for diabetes mellitus and the associated diabetic foot syndrome
- Front reinforcement shell for stable hold of the foot
- Optimum wearing comfort thanks to two individual inflatable air pads in the inner shoe
- Basic body made of durable, light plastic
- Comfortable wear, thanks to the foam padding
- Correct fixing thanks to the circular Velcro® fastening
- Optimized tilting sole
- Optional padding for further adaptation to the leg

REF.530

- Adaptable foam inserts for partial pressure relief from the sole
- Special stockings for diabetes mellitus and the associated diabetic foot syndrome
- Reinforcing shells on the front and back, for increased stability
- Inflatable air pads for optimum fit to the limb and improved stability
- Structure in lightweight, sturdy plasti

ADAPTING FOR THE DOCTOR / TECHNICIAN REF.520/530

- Open all the straps and spread them apart (on the back). Fasten them temporarily on themselves. Remove the padding inside the plastic boot and open it.
- Put on the padding so that the heel rests firmly up against the rear section and then close using the Velcro® fastening (Fig. A).
- Put the plastic boot on the leg from the bottom, spreading the plastic splint (Fig. B).
- Position it so that the splint is vertical and centered on the ankle (Fig. C).
- Close the Velcro® straps, starting from those closest to the toes (Fig.D) and moving upward (Fig. E).
- Inflate the air pads (Fig. F and G) and ensure comfortable fit.
- To inflate: Insert the small pump (Fig. F) with the opening IN into the desired valve and inject air by using the pump. To deflate: Insert the small pump with the opening OUT into the desired valve. It is then possible to use the pump actively to pump the air out or leave it in place and press on the padding until it is completely deflated. The healthcare personnel assigned should teach the patient how to obtain adherence, comfort and adequate compression, for when the boot is reapplied later by the patient.
- If necessary, use the additional padding supplied to adapt the splint to the patient's leg (Fig. H).

Important: Remove the boot only when prescribed by the healthcare operator (to remove the boot follow the above sequence in reverse).

ADAPTING FOR THE DOCTOR / TECHNICIAN REF.505 (Inserting the insert)

- Remove the inner padding and the place the insert into the inner sole of the walker for fit without removing the adhesive protection foil (Fig.1).
- Adapt the insert to the patient as required.
- Remove the adhesive protective foil (Fig.2) and stick to fit into the inner sole of the walker and press firmly into place (Fig.3).

PUTTING ON THE APPLIANCE FOR THE PATIENT REF. 505 (Putting on the special stockings)

- Turn the special stocking inside out up to as far as the heel (Fig.4).
- Pull the stocking up to the heel over the foot (Fig.5).
- Pull the special stocking bit for bit over the foot, heel and calf (Fig.6) and smooth until it has no creases (Fig.7).

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