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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Walker with range of motion

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF. 600		
Size	S	M	L
Shoe size	36/39	39/44	44/47
Plantar flexion	0° - 7,5° - 15° - 22,5° - 30° - 37,5° - 45°		
Dorsal flexion	0° - 7,5° - 15° - 22,5° - 30° - 37,5° - 45°		
Colour	black		

fits right and left

MAINTENANCE

- ☒ Do not bleach ☒ No chemical cleaning
- ☒ Do not iron ☒ Do not tumble-dry
- ☒ Washing instructions:
- Padding: Wash by hand in lukewarm water (max. 30°C) with neutral soap; rinse thoroughly. Dry away from heat sources
- Rigid parts: Clean with a sponge soaked in lukewarm water (max. 30°C) and neutral soap. Dry with a cloth.
- Do not dispose of the product or any of its components into the environment..

INDICATIONS

- Malleolar fractures conservative treatment
- Grade 3 ankle sprains
- Post-surgical care for osteosynthesis of malleolar fractures, ankle ligament recontruction and Achilles tendon tenorraphy
- Walking recovery in consolidation delay of tibial fractures
- Charcot foot

CONTRAINDICATIONS

The application of this protector is advised only when it is documented that the fracture is stable and there are acceptable limits of angular and rotational misalignments.

CHARACTERISTICS AND MATERIALS

- Structure made of durable, light plastic
- Non-slip rubber sole
- Inner padding of the sole made of non-toxic Evazote®
- Struts with Velcro for padding grip
- Soft foam pad with double lining in brushed nylon for Velcro® grip
- Velcro® fastening tapes
- Washable padding

ADAPTING FOR THE DOCTOR / TECHNICIAN

- 1 Undo the Velcro® straps and remove the inner padding from the boot. Insert the foot into the padding. Pay attention to ensure that the padding is positioned at the rear of the heel.
- 2 First close the foot section and then close the calf section of the padding using the Velcro® fastener (fig. A). The pad must be snug the full length but the blood circulation may not be hampered.
- 3 Set the joint as described on the next pages.
- 4 Open the shells with both hands and pull the boot over the padding. Align the shells with the centre axis of the calf (fig. B). If need be, additional pads may be used to ensure better hold and optimum comfort for the wearer.
- 5 Secure the pad using the Velcro® faster of the shells. To do this, remove the separating plastic strip by pulling upwards (fig. C).
- 6 Secure the Velcro® straps of the boot. First close the straps at the toe, and then the others one after another, working from bottom to top (fig. D-E).

SETTING THE JOINTS

- 1 Rotate the plastic cover to open (fig. 1).
- 2 Rotate the cover to the Allen head screw and loosen the screw with the appropriate key (fig. 2).
- 3 Setting the dorsal flexion (fig. 3): Rotate the cover until the opening is positioned over the flexion lock and pull out the pin. Rotate the cover and the pin to the opening with the desired degree of flexion (shown on the mechanism) and insert the pin at this position.
- 4 Setting the plantar flexion (fig. 4): Rotate the cover until the opening is positioned over the extension lock and pull out the pin. Rotate the cover and the pin to the opening with the desired degree of extension (shown on the mechanism) and insert the pin at this position.
- 5 Close the plastic cover again.
- 6 If the brace is to be locked fully, set the pins for dorsal and plantar flexion to the same degree and tighten the Allen head screw.

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