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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Brace for foot drop

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient.

To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF.20050				
Size	XS	S	M	L	
Ankle circumference cm	13/17	17/21	21/25	25/29	

Colour: black - fits right and left

SELECTION/SIZE ACCESSORY

Code	REF.20051				
Size	S	M	L		
Midfoot circumference cm	17/21	21/25	25/29		

Colour: black - specify right or left

MAINTENANCE

- ☒ Do not bleach
- ☒ No chemical cleaning
- ☒ Do not iron
- ☒ Do not tumble-dry

Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Flaccid paralysis of the foot
- Light to moderate spastic equinus foot

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- To be worn ONLY with lace-up shoes.
- Low profile.
- Breathable micropore fabric.
- A soft pad attaches the support to the shoe.
- Tensioning elastic to adjust the degrees of flexion.
- Gel padding over the Achilles tendon.
- Fastens with Velcro® straps.
- Additional midfoot band for wearing with bare feet (optional).

PUTTING ON THE APPLIANCE

- Separate the three parts of the support (shoe pad, tensioning elastic and the tibial cuff) and fully open the tibial cuff.
- With the tibial cuff open, rest it above the back of the ankle, with the gel padding over the Achilles tendon and the bend upwards (fig.1a), and wrap it around the ankle using the two Velcro® straps (fig.1b). Leave the ring completely free and uncovered.
- Insert the shoe pad between the laces and the tongue of the shoe (fig.2).
- Put on the shoe and lace it up (fig.3).
- With the tensioning elastic folded in two, pass it through the D-ring from top to bottom (fig.4).
- Attach the folded end of the elastic to the shoe pad (fig.5).
- Finally, attach the elastic to the support and achieve the desired degree of tension by fastening the two Velcro® ends level with the malleoli (fig.6).

AFFIX OPTION: MIDFOOT BAND

CHARACTERISTICS AND MATERIALS

- To be worn with bare feet or open shoes together with the Affix.
- Low profile.
- Soft, breathable micropore fabric.
- Anatomically designed to fit the left or right foot.
- Fastens with Velcro®.

PUTTING ON THE APPLIANCE

- With the tibial cuff open, rest it above the back of the ankle, with the gel padding over the Achilles tendon and the bend upwards (fig.10a), and wrap it around the ankle using the two Velcro® straps (fig.10b). Leave the D-ring completely free and uncovered.
- Fit the band onto the midfoot with the ring on the top (fig.11).
- With the tensioning elastic folded in two, pass it through the D-ring from top to bottom (fig.12).
- Attach the folded end of the elastic to the ring on the midfoot band (fig.13).
- Finally, attach the elastic to the support and achieve the desired degree of tension by fastening the two Velcro® ends level with the malleoli (fig.14).

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