# fisio

REF. 6002/G - REF. 6002/P

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Board for stabilisation and balance

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Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland

Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89

info@orthoservice.com · www.orthoservice.com

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Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13

info@orthoservice.de · www.orthoservice.de

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
## ENGLISH

Please read the instructions carefully and thoroughly


### DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares with sole responsibility that this product is a Class I medical product and has been manufactured in accordance with the basic requirements of the directives 93/42/EEC and 2007/47/EEC. These operating instructions were created in accordance with the specification of the above-mentioned directives. They aim to ensure the safe and correct use of the medical product.

### SELECTION/SIZE

Item	REF. 6002/P	REF. 6002/G
Size	15x31	30x31
Weight	440 g	900 g

Colour: green · universal

The double board can easily be unscrewed to make two individual boards

High mobility in all directions thanks to the rounded base

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### MAINTENANCE

- Clean with a sponge soaked in lukewarm water (30° C) and neutral soap. Dry with a cloth. Dry away from heat sources
- After use dispose the product and any of its components appropriately

### PRECAUTIONS

If in doubt about how to use the product, contact your physician or an orthopaedic technician. The product should not be used in the vicinity of open flames or strong electromagnetic fields.

### WARNINGS

The product may only be used in the event of medical indication or when prescribed by a physician/physiotherapist. Any changes to the structure or adjustment of the product must be prescribed by a physician or a physiotherapist and performed by an orthopaedic technician. Use of the product is recommended for a single patient. With hypersensitive patients, direct contact with skin may cause redness or irritation. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician.

### WARRANTY

This warranty is governed by the statutory provisions of the country in which the product was purchased. Where relevant that country's regulations concerning warranties and the relationship between the seller and buyer shall apply. If you believe that the terms of the warranty may be applicable, first contact the retailer where the product was purchased. Orthoservice will repair or replace all or part of the device and its accessories if defects in materials and workmanship arise during the period of six months from the date of purchase.

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## ENGLISH

### CHARACTERISTICS AND MATERIALS

- Made from stable, break-resistant plastic with a high resilience
- Non-slip surface with nubs
- Rubber base to ensure adhesion to the ground
- Spacer ring (if required) in order to increase the oscillation angle of the board
- Single board: High mobility in all directions thanks to the rounded base
- Dual board: Can be split into two separate boards, if necessary

### INDICATIONS

- Rehabilitation following ankle sprains
- Rehabilitation of ankle instability
- Physical therapy following operations on ankle ligaments and achilles tendon
- Rehabilitation of knee sprains

### CONTRAINDICATIONS

- Currently no known

### ADAPTING / PUTTING ON THE APPLIANCE

#### EXERCISES WITH SINGLE BOARD

The following exercises are examples only; the program of exercises to be performed must be prescribed by the doctor or physiotherapist.

#### Exercise 1.1: Joint stabilisation (fig. A-C)

Position the rehabilitation board so that you can perform bending and stretching movements (forwards and backwards). Place the foot of the injured leg on the board while keeping the torso at a vertical position to the leg. Take half a step forwards or backwards with the upper body kept in an upright position.

#### Exercise 1.2: Proprioception training (fig. D-F)

#### (to be executed with the knee in flexed position of 20°-30°)

Position the large rehabilitation board so that you can perform bending and stretching movements (forwards and backwards). Stand on the board with bare feet at first and try to adopt a correct position. Correct any wrong positions and postures (e.g. fallen arches) by tensing the muscles and pressing the feet against the board to avoid slipping (this is made easier by the ribbing). The heel and big toe must have good contact with the board. Perform the movements by putting some weight on the injured ankle. The board must turn around on its axis, as far as permitted by the ankle's mobility, WITHOUT EXPERIENCING ANY PAIN (the edges do not need to the touch the floor).

#### EXERCISES WITH DUAL BOARD

The following exercises are examples only; the program of exercises to be performed must be prescribed by the doctor or physiotherapist.

#### Exercise 2.1: Strengthening the knee and foot muscles (fig. G-I)

Position the large rehabilitation board so that you can perform bending and stretching movements (forwards and backwards). Stand on the board with bare feet at first and try to adopt a correct position. Correct any wrong positions and postures (e.g. fallen arches) by tensing the muscles and pressing the feet against the board to avoid slipping (this is made easier by the ribbing). The heel and big toe must have good contact with the board. Perform the movements by putting some weight on the injured ankle. The board must turn around on its axis, as far as permitted by the ankle's mobility, WITHOUT EXPERIENCING ANY PAIN (the edges do not need to the touch the floor).

#### Exercise 2.2: Side stabilisation (fig. K-M)

Position the large rehabilitation board so that you can perform inversion and eversion movements (to the side). Stand on the board with bare feet and try to adopt a correct position; the heel and big toe must have good contact with the board. Perform the movements by putting some weight on the injured ankle. The board must turn around on its axis, as far as permitted by the ankle's mobility, WITHOUT EXPERIENCING ANY PAIN (the edges do not need to the touch the floor).

#### Exercise 2.3: Strengthening the ankles (fig. N-P)

Stand on the board with both feet. The toes are positioned along the centre axis of the board. Now raise the heels so that you are standing on your tip toes. Keep your balance.

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