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 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

Board for stabilisation and balance

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product.

We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient.

To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation.


The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation.

In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF. 6002/P	REF. 6002/G
Size	15x31	30x31
Weight	440 g	900 g
Colour	green	

universal

MAINTENANCE

- Clean with a sponge soaked in lukewarm water (30° C) and neutral soap. Dry with a cloth. Dry away from heat sources
- After use dispose the product and any of its components appropriately.

INDICATIONS

- Rehabilitation following ankle sprains
- Rehabilitation of ankle instability
- Physical therapy following operations on ankle ligamentes and achilles tendon
- Rehabilitation of knee sprains

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Made from stable, break-resistant plastic (Polystyrene) with a high resilience
- Non-slip surface with nubs
- Rubber base (PVC) to ensure adhesion to the ground
- Spacer ring (PVC) (if required) in order to increase the oscillation angle of the board
- Single board: High mobility in all directions thanks to the rounded base
- Dual board: Can be split into two separate boards, if necessary

ADAPTING / PUTTING ON THE APPLIANCE

EXERCISES WITH SINGLE BOARD

The following exercises are examples only; the program of exercises to be performed must be prescribed by the doctor or physiotherapist.

Exercise 1.1: Joint stabilisation (fig. A–C)

Position the rehabilitation board so that you can perform bending and stretching movements (forwards and backwards). Place the foot of the injured leg on the board while keeping the torso at a vertical position to the leg. Take half a step forwards or backwards with the upper body kept in an upright position.

Exercise 1.2: Proprioception training (fig. D–F)

(to be executed with the knee in flexed position of 20°–30°)

Position the large rehabilitation board so that you can perform bending and stretching movements (forwards and backwards). Stand on the board with bare feet at first and try to adopt a correct position. Correct any wrong positions and postures (e.g. fallen arches) by tensing the muscles and pressing the feet against the board to avoid slipping (this is made easier by the ribbing). The heel and big toe must have good contact with the board. Perform the movements by putting some weight on the injured ankle. The board must turn around on its axis, as far as permitted by the ankle's mobility, WITHOUT EXPERIENCING ANY PAIN (the edges do not need to the touch the floor).

EXERCISES WITH DUAL BOARD

The following exercises are examples only; the program of exercises to be performed must be prescribed by the doctor or physiotherapist.

Exercise 2.1: Strengthening the knee and foot muscles (fig. G–I)

Position the large rehabilitation board so that you can perform bending and stretching movements (forwards and backwards). Stand on the board with bare feet at first and try to adopt a correct position. Correct any wrong positions and postures (e.g. fallen arches) by tensing the muscles and pressing the feet against the board to avoid slipping (this is made easier by the ribbing). The heel and big toe must have good contact with the board. Perform the movements by putting some weight on the injured ankle. The board must turn around on its axis, as far as permitted by the ankle's mobility, WITHOUT EXPERIENCING ANY PAIN (the edges do not need to the touch the floor).

Exercise 2.2: Side stabilisation (fig. K–M)

Position the large rehabilitation board so that you can perform inversion and eversion movements (to the side). Stand on the board with bare feet and try to adopt a correct position; the heel and big toe must have good contact with the board. Perform the movements by putting some weight on the injured ankle.

The board must turn around on its axis, as far as permitted by the ankle's mobility, WITHOUT EXPERIENCING ANY PAIN (the edges do not need to the touch the floor).

Exercise 2.3: Strengthening the ankles (fig. N–P)

Stand on the board with both feet. The toes are positioned along the centre axis of the board. Now raise the heels so that you are standing on your tip toes. Keep your balance.

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