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    **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

#### Shock absorber heel cups

##### DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

##### TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

##### SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

##### WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care.

Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician.

The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation.

The orthosis is made and to be used by only one patient.

Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices.In hypersensitive individuals, direct contact with the skin may cause redness and irritation.

In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

##### SELECTION/SIZE

<b>Code</b>	<b>REF. F1350 Talufit 50</b>
<b>Size</b>	<b>2</b> <b>3</b>
<b>Body weight</b>	<b>&lt; 80 kg</b> <b>&gt; 80 kg</b>
<b>Colour:</b>	transparent / Fits right and left

<b>Code</b>	<b>REF. F1352 Talufit 52</b>
<b>Size</b>	<b>1</b> <b>2</b> <b>3</b>
<b>Body weight</b>	<b>&lt; 50 kg</b> <b>&lt; 80 kg</b> <b>&gt; 80 kg</b>
<b>Colour:</b>	anthracite / Fits right and left

##### MAINTENANCE

☒ Do not bleach      ☒ No chemical cleaning  
☒ Do not iron      ☒ Do not tumble-dry  
☒ Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

##### INDICATIONS

- Tallonitis and anterior talalgia
- Sever-Haglund disease
- Achilles tendinopathy
- Calcanear spur
- Retrolcalcaneal Bursitis

##### CONTRAINDICATIONS

Currently no known

##### CHARACTERISTICS AND MATERIALS

- Made of scented and elastic polymer in a double-grid mesh
- Anatomically contoured
- Easily insertable into every shoe
- REF. F1350 only:**

Ideal for those who remain standing or walk for many hours

- REF. F1352 only:**

Ideal for sports activities

##### PUTTING ON THE APPLIANCE

Very easily. Just position them inside your shoe, at the level of the heel. They can enter any shoe, from the most elegant to sportswear.

##### WHAT ARE THEY USED FOR?

They absorb the impact of the foot on the ground, preventing pain and the consequences that these slight but continual shocks may have on both the heel and on the articulations of the leg. Pain in the heel may be harmful to the posture and to the way of walking and these, in their turn may create problems for the entire skeletal structure.

**Ideal for long periods of standing or walking.**

##### WHO IS ADVISED TO USE THEM?

They can be used by anyone. They are also ideal for those who remain standing many hours a day, those who walk a lot and those who practice sport, as prevention for talalgia.

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