
# podifer 20

FLEBOLOGISCHE EINLEGESOHLEN MIT QUERVERSTÄRKUNG - SEMELLE PHLÉBOLOGIQUE AVEC BARRE DE RENFORT - ФЛЕБОЛОГИЧЕСКИЙ СУПИНАТОР - WKŁADKA FLEBOLOGICZNA Z PODPARCIEM POPRZECZNYM - FLEBOLOŠKI ULOŠCI S TRANSVERZALNIM POJACANJEM - PLANTARE FLEBOLOGICO CON BARRA

REF.90020

Take care feel better

REF.90020  
Größe / Size / Taille / Paamep / Rozmiar buta / Velicina / Taglia  
35 36 37 38 39 40 41

FLEBOLOGISCHE EINLEGESOHLEN MIT QUERVERSTÄRKUNG - SEMELLE PHLÉBOLOGIQUE AVEC BARRE DE RENFORT - ФЛЕБОЛОГИЧЕСКИЙ СУПИНАТОР - WKŁADKA FLEBOLOGICZNA Z PODPARCIEM POPRZECZNYM - FLEBOLOŠKI ULOŠCI S TRANSVERZALNIM POJACANJEM - PLANTARE FLEBOLOGICO CON BARRA

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**INDICATIONS:** Vasculopathy · Flatfoot or high arched foot · Metatarsalgia  
These orthopaedic insoles for normal active aid circulation. The special memory foam is very effective and reenergises the capillaries while walking, helping to pump blood back up the legs. The special fabric covering woven with silver thread is antibacterial with odour control.

**ENG. PUTTING ON THE APPLIANCE**  
**FIRST APPLICATION**  
1. Remove the original insole from the patient's shoe.  
2. Place the insole in the shoe, check that it is correctly positioned, and ask the patient to try it on.  
3. Alterations may be made to the insole to suit individual requirements (it is advisable for the patient to wear the insole for a few hours before making any alterations). The front can be shaped with scissors (use the insole that was previously in the patient's shoe for reference). The arched detailing supporting the longitudinal arch can be filed down a little at a time, as can the outer edge of the rear part of the insole.  
4. We recommend removing the insole from the shoe after wear so that it can air.

**SUBSEQUENT WEAR**  
1. Place the insole in the shoe and check that it is correctly positioned.  
2. Place the insole in the shoe and check that it is correctly positioned.  
3. Alterations may be made to the insole to suit individual requirements (it is advisable for the patient to wear the insole for a few hours before making any alterations). The front can be shaped with scissors (use the insole that was previously in the patient's shoe for reference). The arched detailing supporting the longitudinal arch can be filed down a little at a time, as can the outer edge of the rear part of the insole.  
4. We recommend removing the insole from the shoe after wear so that it can air.

**CONTRAINDICATIONS**  
Currently no known contraindications.

**CHARACTERISTICS AND MATERIALS**  
The bottom layer is made of PU and acts as a protective layer for the layers above, providing protection against wear and deformation. It also makes it easier to position the insole in footwear. Lower shell made of rigid PU for the stabilisation and support of the foot in pronation and supination. The presence of the transverse reinforcement supports the metatarsals and cushions the contact with the floor, in addition it strengthens the pump function of the insole, improving the venous return and, at the same time, stimulating the reflex zones for heart and lungs.

**DECLARATION OF CONFORMITY**  
As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). The device has been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

**MAINTENANCE**  
Washing instructions:  
Wash by hand up to 30°C with neutral soap, let it dry away from heat sources. Do not dispose of the product or any of its components into the environment.

**SAFETY PRECAUTIONS**  
It is advisable that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or warts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopaedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

**WARNINGS**  
It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopaedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopaedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be reapplied, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, warts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse reaction, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopaedic product is only guaranteed when all the components are properly in use.