podifer 40

RHEUMA-EINLAGE - SEMELLE RHUMATOÏDE - РЕВМАТИЧЕСКИЙ СУПИНАТОР - WKŁADKA REUMATYCZNA - PLANTARE REUMATICO

REF:90040

Take care feel better

REF:90040

35	36	37	38	39	40	41	42	43	44	45	46
----	----	----	----	----	----	----	----	----	----	----	----

Größe / Size / Taille / Paazep / Rozmiar buta / Velicina / Taglia

RHEUMA-EINLAGE - SEMELLE RHUMATOÏDE - РЕВМАТИЧЕСКИЙ СУПИНАТОР - WKŁADKA REUMATYCZNA - PLANTARE REUMATICO

RHEUMA-EINLAGE - SEMELLE RHUMATOÏDE - РЕВМАТИЧЕСКИЙ СУПИНАТОР - WKŁADKA REUMATYCZNA - PLANTARE REUMATICO

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OSUK Ltd
Units 47&50 North East BIC
Wearfield
Sunderland Enterprise Park East
Sunderland SR5 2TA
Tel +44 (0) 191 51 66 220
info@orthoservice.uk · www.orthoservice.uk

RO+TEN S.r.l.
Sede legale: Via Marconi De Marchi 7 · I-20121 Milano (MI) · Italia
Sede operativa e amministrativa:
Via Comasina 111 · I-20843 Verano Brianza (MB) · Italia
Tel. 0039 0362 18 51 200
info@roplusten.com · www.roplusten.com

UK Rep: UK RepMed Ltd
Wey Court West Union Road
Farnham GU9 7PT (UK)

MADE IN ITALY · 8DCCT00065INT05.0225

Take care feel better

ENG: PUTTING ON THE APPLANCE

FIRST APPLICATION

1. Remove the original insole from the patient's shoe.
2. Place the insole in the shoe, check that it is correctly positioned.
3. Alterations may be made to the insole to suit individual requirements (it is advisable for the patient to wear the insole for a few hours before making any alterations). The front can be shaped with scissors (use the insole that was previously in the patient's shoe for reference). The arched detailing supporting the longitudinal arch can be filed down a little at a time, as can the outer edge of the rear part of the insole.
4. We recommend removing the insole from the shoe after wear so that it can air.

SUBSEQUENT WEAR

1. Place the insole in the shoe and check that it is correctly positioned.
2. We recommend removing the insole from the shoe after wear so that it can dry.

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Top layer of fabric made with woven silver threads with natural anti-bacterial and odour-inhibiting effect.
- The special EVA layer automatically shapes itself according to the pressure it is exposed to, filling out the areas exposed to lower pressures, thus ensuring increased wearing comfort.
- The middle visco-elastic layer with form memory effect, in an optimised thickness and with damping function.
- The bottom layer is made of PU and acts as a protective layer for the layers above, providing protection against wear and deformation. It also makes it easier to position the insole in footwear.
- Lower shell made of rigid PU for the stabilisation and support of the foot in pronation and supination.

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

MAINTENANCE

Washing instructions:
Wash by hand up to 30°C with neutral soap, let it dry away from heat sources. Do not dispose of the product or any of its components into the environment.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

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