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SINGLE PATIENT
MULTIPLE USE

 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

Cervical traction kit

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Item	REF.3019
Size	Universal with linen bag

Item	REF.3020
Size	Universal with water weight bag

MAINTENANCE

Washing instructions:

- For components head harness (REF. 3021) and linen weight bag (REF. 3023): hand wash at max. 30° C with neutral soap, dry away from heat sources
- Metallic parts: clean with a sponge soaked in lukewarm water (30° C) and neutral soap, dry away from heat sources
- Contact an orthopaedic technician for the replacement of used parts
- Do not dispose of the product or any of its components into the environment

SCOPE OF DELIVERY

SET REF. 3019 CERVITRAC WITH LINEN BAG:

- Head harness for traction (REF. 3021) (fig. G)
- Spreader bar (REF. 3034), pulley, cord and hook (fig. H)
- Linen bag (REF. 3023) (fig. L)

SET REF. 3020 CERVITRAC WITH WATER WEIGHT BAG

- Head harness for traction (REF. 3021) (fig. G)
- Spreader bar (REF. 3034), pulley, cord and hook (fig. H)
- Water weight bag (REF. 3022) (fig. I)

INDICATIONS

- Cervical osteoarthritis
- Straightening of the physiological lordosis cervical
- Cervical discopathies
- Cervicobrachialgia
- Reduction of intervertebral spaces

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Over door clamp and pulley
- Head support structure, universal, white
- Spreader bar, universal, silver
- REF. 3019: Linen bag
- REF. 3020: Water weight bag, universal, Max. weight 9kg
- Nylon cord
- Hook

ADAPTING FOR THE DOCTOR / TECHNICIAN

- Attach the pulley to the top edge of a door (fig. A).
- Pass the rope over the rolls of the pulley (fig. B).
- Tie the rope together (fig. C):
 - with the clamping yoke at the front
 - with the hook at the back
- Have the patient sit with a straight back in a seat close to the door (fig. D).
- Attach the chin support to the face/head, to do so (fig. E):
 - open the Velcro® straps
 - adjust the curved chin rest and the flat head rest
 - fasten the Velcro® straps underneath the chin
 - attach the D rings in the intended strap notches
- Mount the bag with the prescribed weight onto the hooks (fig. F).

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