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 PLEASE READ THESE INSTRUCTIONS CAREFULLY AND SAVE THEM

## SET for vertebral traction

### DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG states under its sole responsibility that this product is a Class I medical device that has been manufactured in accordance with the requirements of Directives 93/42/ EEC and 2007/47/EC. These instructions have been written in conformity with these directives, and their aim is to ensure the safe and proper use of this medical device.

### TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

### PRECAUTIONS




It is advised that the pressure exerted by the product does not act on parts of the body suffering from wounds, swelling or lumps. The product should not be tightened too firmly so as to avoid the development of pressure points or compression of the nerves or blood vessels underneath. If in doubt about how to apply, contact an orthopaedic technician. Please read the composition of the product on the inside label carefully. You should not wear the device in the vicinity of an open flame or strong electromagnetic fields.

### WARNINGS

The product has been designed for the specific indications stated below and must be prescribed by a doctor and fitted by an orthopaedic technician. The orthopaedic technician is responsible for fitting the product and providing information about how to use it safely in accordance with individual requirements. The product must be fitted with the utmost care to ensure its efficacy, tolerability and correct function. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. The manufacturer is no longer liable if the product is used or adapted inappropriately. The orthosis is intended to be used only by a single patient; if it is used for more than one patient, the manufacturer accepts no responsibility or liability in accordance with medical device legislation. Direct contact with the skin may cause reddening or irritation in hypersensitive individuals.

In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician. The orthopaedic efficacy of the product may only be guaranteed when all of its components are used. Maximum capacity of the device 9 kg. ORTHOSERVICE AG will not be held liable for higher weights compared to those indicated.

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Do not dispose of the product or any of its components into the environment.

### INDICATIONS

- Lumbago and lumbosciatica
- Lumbar and lumbo-sacral discopathies
- Lumbar paravertebral contractures
- Moderate disk hernias
- Size reduction of disk spaces

### CONTRAINDICATIONS

No known

### PUTTING ON THE APPLIANCE

- Position the platform besides the bed and adjust the height
- Lay the patient on the bed, on his/her back
- Apply the pelvic strap, checking that the lateral straps are aligned along the lower limbs in case one central weight is used:
  - Introduce the "D" rings of the straps in the arch
  - Tighten the cord to the central hole of the arch
  - Pass the cord above the central washer
  - Tighten the hook to the free edge of the cord and fasten the bag with the relative weight
- In case two lateral weights are used
  - Tighten a cord to each of the two "D" rings
  - Pass each cord above the relative side washer
  - Tighten the hook to the free edge of each cord and fasten the bag with the relative weight

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