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
ORTHOSERVICE
RO+TEN

Take care feel better

EN – The descriptions and images in this document are for example and commercial purposes only. Orthoservice reserves the right to modify them according to its needs.

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OSUK Ltd
Units 47&50 North East BIC
Wearfield
Sunderland Enterprise Park East
Sunderland SR5 2TA
Tel +44 (0) 191 51 66 220
info@orthoservice.uk · www.orthoservice.uk

RO+TEN S.r.l.
Sede legale: Via Marco De Marchi 7 · I-20121 Milano (MI) · Italia
Sede operativa e amministrativa:
Via Comasina 111 · I-20843 Verano Brianza (MB) · Italia
Tel. 0039 0362 18 51 200
info@roplusten.com · www.roplusten.com

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4

5

PLEASE READ THESE INSTRUCTIONS CAREFULLY AND SAVE THEM

Motorised trainer for passive mobilisation

SAFETY PRECAUTIONS

Electrical safety

We strongly recommend that you follow these safety guidelines to avoid any risk of fire, electric shock, burns or injury:

- Use a 230V - 50Hz mains power supply.
- Always unplug the trainer from the power supply when not in use.
- Do not clean the trainer while plugged into the power supply.
- Never leave the trainer unattended while plugged into the power supply.
- Do not attempt to use the trainer if the electric cable or plug appear damaged.
- Never use the trainer if it is wet.
- Do not pull the electric cable to move the trainer.
- Do not use the trainer if sprays or compressed oxygen have been used.
- Always turn off the trainer with the on/off button before unplugging it from the power supply.

WARNING

The trainer is not for outdoor use. Consult your physician before beginning any type of exercise programme. If you experience any joint pain, muscle pain, dizziness, abnormal heart rate, or altered breathing patterns, stop exercising and consult your physician.

The trainer is not designed for use by persons with mental incapacity. Keep children and pets away from the trainer while in use. Children and people with disabilities should only use the trainer under supervision. Do not wear loose-fitting clothes or jewellery (bracelets, ankle bracelets, necklaces) while using the trainer. Ensure that the trainer is always on a safe, stable surface; instability may cause the trainer to tip over or slide which could result in injury.

Do not use the trainer if it is worn, damaged or wet. More specifically, we recommend that you take particular care to ensure the cable and plug are in good condition. Do not use the trainer if it has been subjected to any severe impacts, even if there is no obvious damage.

Do not use the trainer near flammable substances or in areas with a high concentration of oxygen.

Ensure that the trainer is never left plugged into the power supply when not in use.

Unplug the trainer from the power supply before changing the pedals or handles.

Never use the trainer under blankets, cushions, mattresses, or layers of paper or cardboard: overheating could cause a short circuit, fire, burns or injury.

Do not use the trainer if the ventilation holes at the sides are blocked. If the trainer does not function normally, unplug from the power supply, without opening or tampering with it, and contact your dealer. Repairs may be made by the authorised service centres recommended by the dealer. Keep this instruction manual in a safe place for future reference as and when needed.

TECHNICAL FEATURES

Power: 230V-50Hz

Electric motor: 230V-50Hz

RPM: 48 min

Motor power: 50W

Weight: 5,7kg

Dimensions: 440x425x285mm

Self-locking device in case of motor overload.

Base supports with holes for using the trainer in the vertical position.

Top handle.

The polarized electrical plug reduces the risk of electric shocks and/or accidents.

Lateral stabilisers with nonslip feet keep the trainer securely on the work surface.

PACKAGING

A Two pedals*

B Ventilation holes

C Electronic control with cable

D On/off button

E Two base supports*

F Two handles*

* These pieces require assembly. The packaging includes a key for assembly.

THE SYSTEM

SPASM DETECTION SYSTEM

The trainer is fitted with a "spasm detection system" to ensure you can exercise safely and effectively. A spasm may occur during exercise which suddenly contracts your muscles.

The Bi-bike recognises the spasm and stops turning for a few seconds, then changes direction to release the spasm. If a second spasm occurs, the trainer locks and emits a beep.

To begin exercising again, just restart the Bi-bike with the start/pause button.

OVERLOAD DETECTION SYSTEM

The trainer is fitted with a safety device that automatically stops the motor it is overloaded or used incorrectly. A beep is emitted to warn the user.

Just press the on/off button to start using the trainer again. The beeping stops automatically after one minute regardless of whether the trainer is started again.

INSTRUCTIONS

- 1 Choose where you wish to use the trainer. For legs: sit on a sofa or chair, or lie down on a bed or on the floor in the gym, with the trainer resting on the floor or on a stable surface. For arms: place the trainer on a table or work surface.
- 2 Choose whether you wish to use the pedals or the handles (fig. 1).
- 3 Slide the pin of the pedal or the handle into the hole and secure by screwing the bolt.
- 4 Check that the power cable is not too taut or tangled.
- 5 Connect the plug to a 230/50 power socket.
- 6 Switch on the trainer with the on/off button (when the button is in the "on" position it lights up to indicate that the trainer is operational).
- 7 Select a type of exercise from the display on the electronic control.
- 8 Press the start/pause button to start, pause or stop exercising.
- 9 At the end of the session, switch off the trainer using the on/off button.

The trainer is programmed for sessions that last a maximum of 15 minutes, which can also be programmed in segments of 3-6-9-12 minutes either on the trainer or using the electronic control. If you wish to exercise for longer, the trainer needs to be switched back on again after it has automatically switched off. While exercising (fig. 2), we recommend that you adopt the correct posture to improve blood flow.

FEATURES AND MATERIALS

- Spasm detection with automatic locking system and warning beep.
- Electronic control with on-board computer.
- 4 adjustable levels of pedal speed:
 - 15 RPM (low) – 30 RPM (medium-low)
 - 45 RPM (medium-high) – 60 RPM (high)
- 5 adjustable levels of duration: 3-6-9-12-15 minutes.
- 2 selectable types of direction: forward or reverse with a two-second pause for a change of direction.
- The trainer has easy-to-attach handles and pedals.
- Rev counter.
- Calorie counter.
- Countdown.
- Preprogrammed exercise function.
- The electronic display shows data and functions: set time, RPM, actual exercise time, number of revolutions completed, calorie consumption, and activation of the preprogrammed exercise function.

INDICATIONS

- General rehabilitation.
- Exercise for elderly people.
- Helps to support the venous pump and lymphatic drainage.
- Complementary exercise for treatments with traditional kinesiotherapy in serious cases of chronic and progressive disability.

CONTROL WITH CABLE

Start/Pause

The start, pause and stop button.

Display

Timer

Use this button to select the duration of your exercise session in minutes (3-6-9-12-15). The trainer automatically stops at the end of the programmed time.

Speed

This button adjusts the RPM speed using four different levels (15 RPM low, 30 RPM medium-low, 45 RPM medium-high, and 60 RPM high).

Auto

Use this button to select preprogrammed exercise. Forward/reverse: the button sets the direction in which the pedals turn (forwards or backwards). The trainer pauses for two seconds before changing direction.

Mode

This button enables you to view data about calorie consumption, total exercise time and the number of RPM on the display.

AUTO – Preprogrammed exercise. Press the AUTO button to select and deselect the preprogrammed exercise function.

If the symbol is underlined, preprogrammed exercise has been selected.

READING THE DISPLAY

AUTO · Preprogrammed exercise.

Press the AUTO button to select and deselect the preprogrammed exercise function.

If the symbol on the display is underlined, preprogrammed exercise has been selected.

15 30 45 60 RPM · Revolutions per minute.

The underlined number on the display corresponds to your selected speed.

3 6 9 12 15 MIN · Session length (in minutes).

The underlined number on the display corresponds to your selected duration.

PLEASE NOTE: all settings may be changed while you are exercising.

0 COUNT (fig. 3) · Revolutions made by the trainer.

Press the MODE button to move to the next function.

0.0 CAL (fig. 4) · Calculate how many calories are burned while exercising.

Press the MODE button to move to the next function.

TIME (fig. 5) · Displays the active exercise time on the trainer.

This can be changed by pressing the TIMER button until the line reaches your chosen time.

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