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RADIO LUCENT

PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

**Immobilizer for arm and shoulder, with elbow support (pediatric)**

**DECLARATION OF CONFORMITY**

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

**TRADEMARKS OF MATERIALS**

Velcro® is a registered trademark of Velcro Industries B.V.

**SAFETY PRECAUTIONS**

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

**WARNINGS**

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

**SELECTION/SIZE**

Code	REF.OFT305
Size	universal
Length of the bag in cm	23 cm
Height	22 cm
Colour	black

fits right and left

**MAINTENANCE**

- Do not bleach  No chemical cleaning
- Do not iron  Do not tumble-dry
- Washing instructions: hand wash at max. 30°C with neutral soap. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

**INDICATIONS**

- Compound humerus fracture
- Glenohumeral subluxation
- Shoulder sprain
- Post-surgical care of humeral fractures and interventions on humeral subluxations
- Postoperative care of the long head of the biceps tendon

**CONTRAINDICATIONS**

Currently no known

**CHARACTERISTICS AND MATERIALS**

- Cotton blend fabric sling
- Soft padded fabric immobiliser swathe
- Shoulder straps with Velcro® fastenings
- Elbow support

**PUTTING ON THE APPLIANCE**

- 1 Insert the arm of the shoulder to be immobilised into the arm pocket.
- 2 Place the straps over the shoulders, cross them at the back, then guide them to the front across the arm pocket and close them with the Velcro® closure. The forearm has to be horizontal, parallel to the ground.
- 3 Attach the fixation band using the corresponding Velcro® closure at the arm pocket, guide around the chest (across the arm to be immobilised and under the opposite arm) and close the band again with the Velcro® closure.

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Velcro® è un marchio di fabbrica registrato di Velcro Industries B.V.