 LŐDZU, UZMANÍGI IZLASIET UN SAGLABÁJIET INSTRUKCIJAS

Rokas un pleca imobilizators ar elkoņa kabatu (bērnū)

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 **DESIGNED AND DEVELOPED IN SWITZERLAND**

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EN – The descriptions and images in this document are for illustrative and commercial purposes only. Orthoservice reserves the right to make modifications without prior notice based on their needs.

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OSUK Ltd
Units 47&50 North East BIC
Wearfield
Sunderland Enterprise Park East
Sunderland SR5 2TA
Tel +44 (0) 191 51 66 220
info@orthoservice.uk · www.orthoservice.uk

RO+TEN S.r.l.
Sede legale: Via Marco De Marchi 7 · I-20121 Milano (MI) · Italia
Sede operativa e amministrativa:
Via Comasina 111 · I-20843 Verona Brianza (MB) · Italia
Tel. 0039 0362 18 51 200
info@roplusten.com · www.roplusten.com

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 PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Immobilizer for arm and shoulder, with elbow support (pediatric)

DECLARATION OF CONFORMITY
As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medcal device.

TRADEMARKS OF MATERIALS
Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS
We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames or strong electromagneticfields. Do not apply in direct contact with open wounds.

WARNINGS
It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient.

To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices.In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

Code	REF.OFT305
Size	universal
Length of the bag in cm	23 cm
Height	22 cm
Colour	black
fits right and left	

MAINTENANCE
⚠ Do not bleach ☒ No chemical cleaning
🧺 Do not iron ☒ Do not tumble-dry
🧼 Washing instructions: hand wash at max. 30°C with neutral soap. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Compound humerus fracture
- Glenohumeral subluxation
- Shoulder sprain
- Post-surgical care of humeral fractures and interventions on humeral subluxations
- Postoperative care of the long head of the biceps tendon

CONTRAINDICATIONS
Currently no known

CHARACTERISTICS AND MATERIALS

- Cotton blend fabric sling
- Soft padded fabric immobiliser swathe
- Shoulder straps with Velcro® fastenings
- Elbow support

PUTTING ON THE APPLIANCE

- Insert the arm of the shoulder to be immobilised into the arm pocket.
- Place the straps over the shoulders, cross them at the back, then guide them to the front across the arm pocket and close them with the Velcro® closure. The forearm has to be horizontal, parallel to the ground.
- Attach the fixation band using the corresponding Velcro® closure at the arm pocket, guide around the chest (across the arm to be immobilised and under the opposite arm) and close the band again with the Velcro® closure.

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