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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Wrist immobilizer (pediatric)

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medcal device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagneticfields. Do not apply in direct contact with open wounds.

WARNINGS

It is advised that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient.

To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation.

The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices.In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF.OFT500
Size	universal
Length cm	15,5
Color	blue

specify right or left

MAINTENANCE

 Do not bleach No chemical cleaning
 Do not iron Do not tumble-dry
 Washing instructions: hand wash at max. 30°C with neutral soap. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Wrist distortion
- Compound fracture of the radius and ulna extremities
- Wrist tendonitis
- Consequences of fracture of the wrist or after early removal of the plaster

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Aluminium palm splint, preformed and malleable (removable)
- Double-lined foam design
- Reinforcing dorsal stays
- Fastening with Velcro® straps and buckles

PUTTING ON THE APPLIANCE

- Open all the Velcro® straps and pull the straps out of the buckle (fastening them provisionally on themselves).
- Position the wrist so that the metal stay is on the same side as the palm of the hand.
- Fixing the brace:
- Pass the straps through the rings provided, tighten and fasten with Velcro®: the central one and the one for the forearm.
- Fasten the Velcro® hand strap.

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