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UK **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM****Wrist immobilizer with thumb grip (pediatric)****DECLARATION OF CONFORMITY**

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF.OFT502
Size	universal
Length cm	16
Colour	blue

specify right or left

MAINTENANCE

- Do not bleach No chemical cleaning
- Do not iron Do not tumble-dry
- Washing instructions: hand wash at max. 30°C with neutral soap. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Wrist sprain
- Wrist tendonitis
- Consequences of wrist fracture or after early removal of the plaster
- Fastening with Velcro® straps and buckles

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Aluminum thumb splint preformed and malleable (removable)
- Aluminum palm splint preformed and malleable (removable)
- Double-lined foam padding
- Reinforcing dorsal plastic stays
- Fastening with buckles and Velcro® straps

PUTTING ON THE APPLIANCE

- 1 Open all the Velcro® straps and pull the straps out of the buckle (fastening them provisionally on themselves)
- 2 Apply the brace such that the metal splint is positioned at the palm and the thumb is inserted into the provided support
- 3 Fixing the brace:
 - Wrap the central strap around the wrist, pass it through the relevant buckle, tighten it and fasten it with the Velcro® strap.
 - Pass the remaining straps through the relevant buckles, tighten them and fasten them with the Velcro®, first the hand one and then the forearm one.
- 4 After fixing the brace fasten the Velcro® strap of the thumb opening

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