ORTHO14-100

ORTHO14-105

ORTHO14-310

ORTHO14-305

ORTHO14-315

ORTHO14-500

ORTHO14-502

ORTHO14-700/2
ORTHO14-700/3

ORTHO14-805

ORTHO14-800

ORTHO14-810

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LATEX FREE

RADIO
LUCENT

PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Bivalve ankle brace (pediatric)

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient.

To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation.

The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF.OFT800
Size	universal
High cm	16
Colour	black

fits right and left

Maintenance

- Do not bleach
- No chemical cleaning
- Do not iron
- Do not tumble-dry
- Washing instructions: hand wash at max. 30°C with neutral soap. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Grade II and III sprains
- Functional recovery phase after surgical treatment
- Functional recovery after malleolar fractures
- Chronic instability of TT and sporting activity

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Double-lined foam padding
- Anatomically shaped lightweight durable plastic bivalve structure
- Velcro® adjustment fits any heel width
- Easy to wear thanks to full opening

PUTTING ON THE APPLIANCE

- 1 Open all Velcro® straps.
- 2 Ensure that the inner inflatable air pads are positioned correctly on the plastic structure.
- 3 Position the heel element under the heel in such a way that the rounded section of the heel element is pointing to the front of the foot.
- 4 Position the support shells at the sides of the ankle bone; if necessary adjust the width of the heel element using the Velcro® closure designed for this purpose, and adjust the width of the heel element.
- 5 Close the orthosis by pulling firmly on the Velcro® strap.

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Velcro® est une marque enregistrée de Velcro Industries B.V.