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REF.OFT805

Walker (pediatric)

DESIGNED AND DEVELOPED IN SWITZERLAND

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EN – The descriptions and images in this document are for illustrative and commercial purposes only. Orthoservice reserves the right to make modifications without prior notice based on their needs.

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OSUK Ltd
Units 47&50 North East BIC
Wearfield
Sunderland Enterprise Park East
Sunderland SR5 2TA
Tel +44 (0) 191 51 66 220
info@orthoservice.uk · www.orthoservice.uk

RO+TEN S.r.l.
Sede legale: Via Marco De Marchi 7 · I-20121 Milano (MI) · Italia
Sede operativa e amministrativa:
Via Comasina 111 · I-20843 Verano Brianza (MB) · Italia
Tel. 0039 0362 18 51 200
info@roplusten.com · www.roplusten.com

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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Walker (pediatric)

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medcal device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames or strong electromagneticfields.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices.In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF.OFT805		
Size	M	L	XL
Shoe size	23/25	26/28	29/31
Height from ground cm	22	25	27
Colour	black		
fits right and left			

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Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Conservative treatment of malleolar fractures
- Severe sprains (3rd degree)
- Post-surgical phase after osteosynthesis of ankle fractures, stabilisation of ligament injuries of the ankle and Achilles tendoni
- Stable foot, ankle, tibia and distal fibula fractures

CONTRAINDICATIONS

The application of this protector is advised only when it is documented that the fracture is stable and there are acceptable limits of angular and rotational misalignments.

CHARACTERISTICS AND MATERIALS

- Structure made from lightweight, durable plastic
- Non-slip rubber sole
- Inner padding of the sole is made from non-toxic Evazote®
- Struts held in place with Velcro® for a firm grip on the padding
- Soft foam pad with double lining in brushed nylon for Velcro® fastening
- Velcro® fastening straps

PUTTING ON THE APPLIANCE

- Undo all the straps and fasten them back on themselves for the moment
- Remove the padding from the boot. Slide your foot into the padding so that the back fits snugly around your heel.
- Fasten the padding and Velcro® over the foot and then around the leg (fig. A). The entire length of the padding should fit snugly, but not so tight that it restricts circulation.
- Open the struts of the boot with both hands and slip the boot on over the padding, aligning the struts with the central line of the ankle
- Extra padding can be added if necessary to ensure a better and more comfortable fit.
- Fasten the boot with the straps; start with the strap over your foot and work your way up the leg

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