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ORTHO14-810

REF.OFT810

Shower for lower leg and foot (pediatric)

DESIGNED AND DEVELOPED IN SWITZERLAND

ORTHOSERVICE
RO+TEN
Take care feel better

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Distributed by
Headquarter: ORTHOSERVICE AG
Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland
Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89
info@orthoservice.com · www.orthoservice.com

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OSUK Ltd
Units 47&50 North East BIC
Wearfield
Sunderland Enterprise Park East
Sunderland SR5 2TA
Tel +44 (0) 191 51 66 220
info@orthoservice.uk · www.orthoservice.uk

RO+TEN S.r.l.
Sede legale: Via Marco De Marchi 7 · I-20121 Milano (MI) · Italia
Sede operativa e amministrativa:
Via Comasina 111 · I-20843 Verano Brianza (MB) · Italia
Tel. 0039 0362 18 51 200
info@roplusten.com · www.roplusten.com

SINGLE PATIENT
MULTIPLE USE

RADIO
LUDENT

LOW
PROFILE

ORTHOSERVICE
RO+TEN
Take care feel better

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UK **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM****Shower for lower leg and foot (pediatric)****DECLARATION OF CONFORMITY**

ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF.OFT810				
Size	1.	2.	3.	4.	5.
Calf circ. cm	15/19	19/23	22/24	23/27	25/29
Length of foot cm	10	11	12	13	15
Splint height cm	16	19	23	25	28
Farbe	semitransparent				

specify right or left

MAINTENANCE

- Do not bleach No chemical cleaning
- Do not iron Do not tumble-dry
- Washing instructions: Hand wash up to 30°C with mild soap (using a sponge is recommended); leave to dry away from direct heat (or better still dry with a cloth).

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Foot drop
- Mild to moderate spastic equinus

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Polypropylene splint
- Can be heated and moulded to fit
- Fastens with a Velcro® strap around the calf

PUTTING ON THE APPLIANCE

- 1 If necessary, trim the splint and mould to fit with the help of a heat gun.
- 2 Don the splint and then tighten and fasten the Velcro® strap around the leg.

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Velcro® è un marchio di fabbrica registrato di Velcro Industries B.V.