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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Shower for lower leg and foot (pediatric)

DECLARATION OF CONFORMITY

ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF.OFT810				
Size	1.	2.	3.	4.	5.
Calf circ. cm	15/19	19/23	22/24	23/27	25/29
Length of foot cm	10	11	12	13	15
Splint height cm	16	19	23	25	28
Farbe	semitransparent				
specify right or left					

MAINTENANCE

- Do not bleach
- Do not iron
- No chemical cleaning
- Do not tumble-dry
- Washing instructions: Hand wash up to 30°C with mild soap (using a sponge is recommended); leave to dry away from direct heat (or better still dry with a cloth).

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Foot drop
- Mild to moderate spastic equinus

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Polypropylene splint
- Can be heated and moulded to fit
- Fastens with a Velcro® strap around the calf

PUTTING ON THE APPLIANCE

- If necessary, trim the splint and mould to fit with the help of a heat gun.
- Don the splint and then tighten and fasten the Velcro® strap around the leg.

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