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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Padded dorso-lumbar orthopedic orthosis with a reinforced structure

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.
Sensitive® is a registered trademark of EUROJERSEY

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Item	REF. 33009					
Size	XS	S	M	L	XL	XXL
Pelvic circ. cm	70/80	80/90	90/100	100/110	110/120	120/130
Height at front cm	16	16	16	16	16	16
Back height cm	50	50	50	52	52	52
Colour	grey					

OPTION: frame REF.33009/T, size 1-L. This special frame fits sizes XS to M - with XXL back height

MAINTENANCE

- ☒ Do not bleach
- ☒ No chemical cleaning
- ☒ Do not iron
- ☒ No dryer or tumble-dry
- ☒ Washing instructions:

- Slide out the metal part
- Hand wash in lukewarm water with mild soap; rinse thoroughly.
- Leave to dry away from direct heat.
- Slide the metal part back into place

Dispose of this product and its parts responsibly

INDICATIONS

- Back pain caused by arthrosis or kyphosis.
- Recovery from fractures of the thoracic and/or upper lumbar vertebrae.
- Support for thoracic vertebral compression fracture caused by osteoporosis or metastasis.
- Other conditions as prescribed by the doctor.

KONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Made of reinforced CarbonTrace, a stretchy, breathable, ladder proof, antistatic and bacteriostatic fabric reinforced with carbon fibre
- Adjustable padded shoulder straps with non-release Velcro® fastening system
- Imbottiture aggiuntive per le spalle
- The elasticated rear panel of the shoulder straps can be positioned at any height
- Adjustable padded shoulder straps with non-release Velcro® fastening system
- Front Velcro® fastening system
- The abdominal section is reinforced with angled spiral stays
- Mouldable inner frame in aluminium alloy

PUTTING ON THE APPLIANCE

DONNING THE ORTHOSIS FOR THE FIRST TIME

- The orthosis comes with a rigid frame already positioned in the pocket; if this has not been prescribed, open the bottom of the pocket, slide out the frame (fig. A) and reclose the pocket with the Velcro®. If required, shape the frame according to the patient's anatomy.
- Undo the orthosis and attach the Velcro® straps to the side strips for the time being (fig. B). The ends of the shoulder straps can also be attached to the side strips for the moment (fig. B) or folded back on themselves.
- Put on the orthosis by sliding the arms through the shoulder straps as if it were a backpack (fig. C).
- Close the orthosis at the front securely with the central Velcro® fastening (fig. D). Check that the front fastening is properly centred.
- Start to tighten the orthosis:
 - Attach the lower straps to the Velcro® strip at the front (fig. E).
- Remove the elasticated panel from the back and reattach with Velcro® at the desired height (fig. F).
- Adjust the shoulder straps:
 - Position the ends on the Velcro® strip at the front and undo the straps from the corresponding Y-Velcro® (fig. G) fasteners.
 - Grasp the straps and pull towards the front (fig. I).
 - Cut off any excess and attach them to the corresponding Y-Velcro® fasteners (fig. M).
- If necessary, slide the additional pads onto the shoulder straps until they are in the required position

SUBSEQUENT WEAR

- Undo the orthosis and attach the Velcro® straps to the side strips for the time being (fig. N). The ends of the shoulder straps can also be attached to the side strips for the moment (fig. N) or folded back on themselves.
- Put on the orthosis by sliding the arms through the shoulder straps as if it were a backpack (fig. O).
- Close the orthosis at the front securely with the central Velcro® fastening (fig. P). Check that the front fastening is properly centred.
- Start to tighten the orthosis:
 - Attach the lower straps to the Velcro® strip at the front (fig. Q).
- Grasp the straps and pull towards the front (fig. R). Attach them to the Velcro® on the strip at the front (fig. S).

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