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
SINGLE PATIENT
MULTIPLE USE

 **PLEASE READ THE INSTRUCTIONS CAREFULLY AND THOROUGHLY AND KEEP THEM FOR FUTURE REFERENCE**

Medical compression stockings

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

MAINTENANCE

-  Do not bleach  Do not dry clean
-  Do not iron  Do not tumble dry
-  Washing instructions: Hand wash up to 30°C with mild soap; rinse without wringing; dry away from direct heat.

If the stocking is looked after properly, the medical efficacy of Orthoservice therapeutic stocking is ensured for six months.

IMPORTANT

- Grip-top design: lotions, talc, sweat or normal body oils can affect the silicon top bands on stockings and lessen their grip. Clean with water up to 30°C and soap as and when necessary.
- One-leg tights: ensure the Velcro® is properly fastened (so the hoop and loop tapes are fully engaged) before washing to prevent the tights from snagging or laddering or damage to the Velcro®.

Dispose of this product and its parts responsibly.

PRECAUTIONS

Do not wear near fire or flames. Not recommended for wearing while sleeping at night. Compression stocking can be worn only after a wound has been treated. Compression stocking must never be worn with and/or over an open wound.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

INDICATIONS

(THE COMPRESSION CLASS IS AT THE DISCRETION OF A DOCTOR)

- Functional disorders, such as heavy, tired legs, without visible varicose veins
- Mild varicose veins associated with subjective disorders (spider veins, reticular veins)
- Varicose veins during pregnancy with or without oedema
- Localised oedema in the evening (venous)
 - Primary varicosities with oedema
 - Advanced primary varicosities with oedema (main vein and branch varicosities)
- Chronic venous insufficiency (CVI)
- Prophylaxis of thrombosis and embolisms
- Thrombophlebitis or venous thrombosis (superficial thrombosis) possibly with oedema
- DVT (deep vein thrombosis)
- Post-thrombotic oedema with or without secondary follow-up therapy
- To ensure successful therapy:
 - following sclerotherapy
 - following an elective outpatient procedure for various veins
 - following stripping
 - Before sclerotherapy or surgery to prevent the onset of oedema
- Ulcer prevention and therapy
 - après un stripping
 - Postoperative oedema
 - Post-traumatic oedema after basic therapy
- Reversible lymphoedema following complex decongestive therapy
- Irreversible lymphoedema with skin becoming hard and tight following complex decongestive therapy
- Lipoedema of grade II and above
- Idiopathic cyclic oedema
- Lymphoedema following basic therapy
- Congestive disorders caused by immobility
- Angiodysplasia

CONTRAINDICATIONS

- Advanced peripheral arterial occlusive disease
- Decompensated heart failure
- Septic phlebitis (venous inflammation caused by microorganisms)
- Phlegmasia coerulea dolens
- Weeping dermatosis
- Intolerance to the material used for graduated compression stocking
- Loss of sensation in the extremities
- Severe peripheral neuropathy (e.g. diabetes mellitus)
- Chronic polyarthritits

The company accepts no responsibility for improper use of this product.

FEATURES AND MATERIALS

- Tubular stretchy medical compression stocking
- Made on new-generation computerised circular knitting machines
- Any seams are flat and extremely fine
- Made with the finest quality Swiss yarns:
 - fine
 - soft
 - guaranteed to be long-lasting and resistant to ointments and creams
 - excellent for maintaining compression and elasticity over time
 - Yarn-dyed for more even colours and better long-term elasticity
 - A special weave is used to make the toe and heel even stronger
 - Extensive quality control involves visual inspection and testing by equipment with sensors
- Unisex

QUALITY CONTROL

This product has undergone meticulous quality control testing before being packaged. If you notice a clear manufacturing fault, please return the product complete with its packaging to the store where you purchased it.

PUTTING ON THE APPLIANCE

It is advisable to put on this medical stocking before getting out of bed in the morning; legs must be dry and preferably sprinkled with a little talc. Wearing rubber gloves can help to slide the stocking more easily over the legs and keep the tension evenly distributed thereby ensuring they fit correctly. If there is a dressing on the leg, it is advisable to first put on a small pair of nylon tights as this will help the medial stocking to slide on more easily.

- A. [Open-toe models only] Simply place the foot slip over the foot, and check it is properly positioned (fig.1).
- B. Orthoservice socks and tights can be put on just like normal socks and tights (fig.2). They do not need to be rolled or turned over beforehand but can be simply bunched up as normal.
- C. Put on the stocking evenly, keeping your hands flat and without forcing the material, up the length of the leg to be covered (fig.3-4): ankle, calf, thigh and groin (depending on the model). Avoid creating any localised tension and the medical stocking will provide the perfect level of compression. The toes will be free if wearing an open-toe model. Take care to always position the heel correctly.
- D. [One-leg-tights only] Fasten the Velcro® at the waist.
- E. [Open-toe models only] Remove the foot slip. (fig.5).
- F. How to correctly position Orthoservice compression stockings and panty stockings:
 - a. Calf-length stockings: the top edge should be positioned below the knee.
 - b. Mid-thigh stockings: the compression fabric covers half the thigh.
 - c. Panty stockings and tights: the compression fabric should come up to the groin.

HOW TO REMOVE YOUR STOCKING

- A. Rubber gloves can make it easier to remove your stocking.
- B. Gently grip the top of the stocking and slide downwards.
- C. When you reach the heel, slide in your thumbs, widen the fabric away from the foot and remove.

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