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 PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Dynamic bunion splint

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medcal device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames or strong electromagneticfields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices.In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF. 21361
Size	Universale
Colour	white
Fits right and left	

Code	REF. 21361
Size	Universale
Colour	white
Fits right and left	
MAINTENANCE	
⚠ Do not bleach	⚠ No chemical cleaning
⚠ Do not iron	⚠ Do not tumble-dry
🧺 Washing instructions: Hand wash in lukewarm water with mild soap; rinse thoroughly. Leave to dry away from direct heat.	
Do not dispose of the product or any of its components into the environment.	

INDICATIONS

- Hallux valgus (bunion).
- Mild or moderate malpositioning of the big toe
- Bursitis in the metatarsophalangeal joint of the big toe

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- The polymer splint is hinged at the metatarsophalangeal joint.
- Soft foam padding.
- Adjustable Velcro® straps.
- Suitable for wearing at night, barefoot, or in comfortable shoes

PUTTING ON THE APPLIANCE

- Open out the straps to their full extent but without sliding them out of their loops.
- Slide the splint over the foot so that the plastic part is aligned with the side of the big toe and the hinge is lined up with the metatarsophalangeal joint (Fig. A).
- Tighten the foot strap and fasten with Velcro® (Fig. B).
- Tighten the big toe strap and fasten with Velcro® (Fig. C).

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