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 PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Carbon-fibre ankle foot orthosis

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medcal device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF.KTAFO		
Size	S	M	L
Height cm	36	39	42
Shoe size	33-38	38-43	43-48
Colour	black		
specify right or left			

MAINTENANCE

-  Do not bleach
-  Do not iron
-  No chemical cleaning
-  No tumble-dry
-  Washing instructions: Padding- Wash by hand in warm water (30°) with neutral soap; rinse properly. Let it dry away from sources of heat. Rigid part (in carbon fibre) - Rub with a sponge soaked in warm water and neutral soap. Dry with a cloth.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Limping foot
- Ictus
- Multiple sclerosis
- Neurological illnesses
- Amputation of forefoot

CONTRAINDICATIONS

- Macroscopic deformity of the foot and ankle
- Acute edema
- From moderate to severe spasticity
- Contracture of front and rear tibial muscles, that cannot be healed with the simple apparatus posture

CHARACTERISTICS AND MATERIALS

- Spiral structure for resistance and elasticity; pre-loaded to provide constant support to the foot
- Structure in reinforced carbon fibre for excellent support. Sole in thin carbon fibre, that can be adapted to different shoes
- Light
- Calf Velcro® closure, with return loop
- Calf padding in MTP finished in technical fabric in contact with the body
- Washable

PUTTING ON THE APPLIANCE

- Open the calf Velcro®, by temporarily folding it on itself.
- Remove the sole inside the shoe (fig.A).
- Introduce the orthosis in the shoe (fig.B).
- Re-introduce the sole inside the shoe above the brace foot support (fig.C).
- Fit the shoe with the inserted orthosis (fig.D).
- Tie the calf Velcro® strap (fig.E+F).

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