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




REF. 21006

## Thumb spica

DESIGNED AND DEVELOPED IN SWITZERLAND

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EN – The descriptions and images in this document are for illustrative and commercial purposes only. Orthoservice reserves the right to make modifications without prior notice based on their needs.

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##### OSUK Ltd

Units 47&50 North East Bic  
Wearfield  
Sunderland Enterprise Park East  
Sunderland SR5 2TA  
Tel +44 (0) 191 51 66 220  
info@orthoservice.uk · www.orthoservice.uk

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 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

## Thumb splica

### DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

### TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

### SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagneticfields. Do not apply in direct contact with open wounds.

### WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient.

To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation.

The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices.In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

### SELECTION/SIZE

Code	REF. 21006	
Taille	S	L
Wrist circumference cm	14-18	18-22
Colour	grey	

Fits right and left

### MAINTENANCE

⚠ Do not bleach ☒ No chemical cleaning  
⚠ Do not iron ☒ Do not tumble-dry  
🗑 Washing instructions: Hand wash in lukewarm water with mild soap; rinse thoroughly. Leave to dry away from direct heat. Dispose of this product and its parts responsibly.

### INDICATIONS

- Rhizarthrosis
- Thumb tendonitis
- De Quervain's tenosynovitis
- Metacarpophalangeal sprain
- A Stener lesion
- The after-effects of Bennet's fracture
- The after-effects of trapezium reconstructive surgery and thumb surgery

### CONTRAINDICATIONS

There are no known contraindications at this time

### FEATURES AND MATERIALS

- A seamless, welded design for maximum comfort
- Anatomically designed
- The aluminium alloy core can be moulded to fit that shape of the patient's hand
- Covered with a strong, padded technical fabric that is breathable and very fine
- Adjustable Velcro® strap for easy and intuitive fastening
- Fits right and left

### DONNING THE SPLINT

- 1 Undo the wrist strap and fasten it back on itself for the time being
- 2 Open out the thumb strap and fasten it with Velcro® (fig. 1)
- 3 Slide the splint onto the hand (fig. 2)
- 4 Check the shape of the splint and adapt it to fit the shape of the hand if necessary
- 5 Slide the wrist strap through the ring provided (fig. 3), tighten it (fig. 4) and fasten with Velcro® (fig. 5)
- 6 Also adjust the strap around the thumb so that it fits snugly

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