


ERSTE ANWENDUNG
FIRST APPLICATION
PREMIÈRE APPLICATION
Первое ПРИЛОЖЕНИЕ

STOSOWANIE PO RAZ PIERWSZY
PUKEMISOJUEET ENSIMAIS-
LES PUJEMISKERALLE
HEY PRODUCT EERSTE GEBRUK
PRIMA ORTOZES ULIKSANA
POSTAVLJANJE ORTOZE
(PRVA PRIMJENA)
PRIMERA APLICACION
PRIMA APPLICAZIONE

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Do not iron

Do not tumble-dry

Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources.

Disposal of this product and its parts responsibly

INDICATIONS

Post-artroscopic care.

Surgical care (compatible with the degree of stabilization required).

Conservative treatment of hip dislocations.

Arthroscopy of the hip joint.

Coxalgia or inoperable joint misalignment.

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

Made of breathable Aer-OS fabric.

Lateral and anterolateral spiral support straps.

Front Velcro® fasteners.

Colour-coded Velcro® areas for fastening the straps.

Elastic calibrated rods, with lateral crossing, for:

Prevention of excessive adduction

Partial limitation of flexion-extension

Prevention of overextension

A rigid bar that locks the joint in flexion and neutral abduction during the initial days following a dislocation

REF.91007N - Adjustable (telescopic) rod for Hipocross 2.0:

limit extension (from 0° to 120° in 10° increments)

limit flexion (from 0° to 120° in 10° increments)

regulate adduction.

PUTTING ON THE APPLIANCE FIRST APPLICATION

1 Position the brace on the hip with the vertical Velcro facing outwards along the median of the body (fig. A)

2 Fasten the Velcro® pelvis strap, checking that the fabric is stretched out properly (fig. B)

3 Fasten the Velcro® thigh strap; hold the loop at the end of the strap if this is easier and check that the fabric is stretched out properly (fig. C)

4 Attach the elastic straps on the lower parts of the pelvis (blue - fig. D) and thigh (orange - fig. E)

5 Cross the elastic straps at the sides (on the trochanter) and fasten: fasten the strap with the orange tip at the front on the pelvis (opposite iliac crest) (fig. F); fasten the strap with the blue tip on the thigh (fig. G). The fit can be tightened or loosened using the hooks at the front and back.

6 Use the bar provided to lock the flexion and abduction of the joint in the neutral position; attach the bar to the side of the Velcro® and fasten the straps provided around the thigh and pelvis. (fig. H)

REF.91007N

If necessary, adjust the height of the adjustable rod (REF.91007N) using the screws provided in the pelvis area and in the thigh area (fig. I).

8 To use the adjustable rod (REF.91007N), attach it to the side strip with Velcro® (fig. L).

9 Fasten the straps around the thigh and pelvis (fig. L-M)

10 The elastic straps can also be fastened over the stay as described above for a tighter fit (fig. O)

ADJUSTMENT OF THE JOINTS

11 Extension adjustment (fig. P)

12 Flexion adjustment (fig. Q)

13 This washer (fig. R) can be turned so that the respective wedge (fig. S) locks in the adjusted flexion-extension setting

14 Hole (fig. T) for the screw provided (to hold the wedge in place)

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SUCCESSIVE APPLICATIONS

1 Position the brace on the hip with the vertical Velcro® facing outwards along the median of the body (fig. 1)

2 Fasten the Velcro® pelvis strap, checking that the fabric is stretched out properly (fig. 2)

3 Fasten the Velcro® thigh strap; hold the loop at the end of the strap if this is easier and check that the fabric is stretched out properly (fig. 3)

4 Hold the elastic straps, cross them at the sides and fasten: fasten the strap with the orange tip at the front on the pelvis (opposite iliac crest) (fig. 4); fasten the strap with the blue tip on the thigh (fig. 5)

5 Apply the brace

PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Brace for hip and coxalgia

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

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Velcro® è un marchio di fabbrica registrato di Velcro Industries B.V.