REF.23070
Lumbobast


REF.23070F
Lumbobast Forte

REF.23070P
• Zusätzliches Pelotten-Kit für Lumbobast • Additional pads kit for Lumbobast • Kit plaques supplémentaires pour Lumbobast
• Дополнительная пластина для Lumbobast • Dodatkowy zestaw podkładek dla Lumbobast
• Kit pelotte agglutivante per Lumbobast

Schnittführung, Cutting guide, Guide de coupe, Guida al taglio, Przewodnik cięcia, направляющая для резки

- Achten Sie darauf, dass der Hakenklettverschluss vollständig am Saum anliegt und bündig mit dem Stoff abschließt
- Attention, veillez à ce que le Velcro® hook completely adheres to the velvet, flush with the fabric
- Attention, veillez à ce que le Velcro® crochet adhère complètement au velours, au ras du tissu
- Attenzione, far aderire completamente il Velcro® uncino al vellutino, a filo del tessuto
- Внимание: крючок-липучка Velcro® должен полностью прилегать к бархату на одном уровне с тканью
- Uwaga, upewnij się, że zaczepiony rzep Velcro® przylega całkowicie do aksamitu, równo z tkaniną

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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

REF.23070 Orthosis with a facilitated traction system to stabilize the lumbar rachis and relieve pressure with rigid thermoformable pads
REF.23070F Orthosis with a facilitated traction system to stabilize the lumbar rachis and relieve pressure with additional pads kit

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Item	REF.23070 - REF.23070F		
Size	1	2	3 (bariatric)
Pelvic circumference	70-106	106-145	130-170
Front height cm	16	16	16
Back height cm	24	24	24

Color: grey

Option: REF.23070P additional pads kit for Lumbofast

MAINTENANCE

- ☒ Do not bleach ☒ No chemical cleaning
 - ☒ Do not iron ☒ Do not tumble-dry
 - ☒ Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources.
- Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Lumbago
- Sciatic lumbalgia and cruralgia
- Spondylarthrosis and lumbar discopathies
- Slight traumas of the lumbosacral column and paravertebral lumbar contractures
- Postoperative recovery
- After-effects of stable compression fractures of the lumbar vertebrae
- Treatment of fracture of transverse processes

ONLY REF.23070F

- Degenerative lumbar scoliosis
- Spondylosis with or without vertebral listhesis

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Bi-coupled base fabric: a three-dimensional mesh on the inside with an exclusive printed technical fabric on the outside. Rigid for maximum hold and breathable for maximum comfort
- Lanyard and pulley traction system for optimum stabilization and maximum traction, with minimum effort
- Lanyard terminals have an anatomical and easy grip
- The rear adjustment system allows you to cut your orthosis to a perfect size
- A rigid thermoformable pad applied with snap buttons provides stability to the column and protects against contact with lanyards.
- Polypropylene reinforcement to protect the lanyard
- Front handles for easy Velcro® closure
- Thermoformable front pads to be inserted in the pockets provided
- Anatomical and modern design
- Lightweight
- Option: thermoformable pads kit (large back and sides)

APPLICATION

FIRST APPLICATION

- Adjust the length of the orthosis:
 - Measure the patient's pelvic circumference
 - Release the buttons and remove the back pad (fig.A1)
 - Remove the back part (fig.A2)
 - Cut off the excess using the pre-cut flattenings and numbers engraved on the orthosis itself as reference (fig.A3). BE CAREFUL to not over-cut (if you are not sure, test first by folding over the excess part and proceed to cut only once verified)
 - Re-attach the back part with Velcro® (fig.A4), re-thread the pad (fig.A5), and re-hook the buttons (fig.A6)
- Open the orthosis (making sure the tie-rods are well loosened/widened) and provisionally fasten the Velcro® ends to the side (fig.B)
- Position the orthosis so that the label is facing upward. Put it on and secure it tightly with the center Velcro® fastener (fig.C). Note: Tighten it well, all the way up.
- Make sure that the Velcro® fastener is well centered.
- Tension the orthosis by pulling the lanyard system through the Velcro® ends (fig.D1) and then hooking them over the front velcraable band (fig.D2).
- If the lanyards are too long, you can shorten them by tying a knot (fig.E1) and heat-cutting the excess part off (fig.E2)
- If necessary, the back pad can be heat-formed according to the patient's anatomy
- The front pads should be inserted on prescription into the appropriate pockets (fig.F).

ONLY FOR Lumbofast Forte REF.23070F

- To apply the large back pad, detach the supplied one, slide it off the orthosis, insert the large one in its place, and fasten it with the snap fasteners (fig.G1 and G2)
- To apply the side pads, slip them in the orthosis and fasten them with Velcro® on the desired spot (fig.H)

SUBSEQUENT APPLICATIONS

- Open the orthosis making sure the tie-rods are well loosened/widened and provisionally fasten the Velcro® ends on the side (fig.I)
- Position the orthosis so that the label is facing upward. Put it on and secure it tightly with the center Velcro® fastener (fig.L). Note: Tighten it well, all the way up
- Make sure that the Velcro® fastener is well centered
- Tension the orthosis by pulling the lanyard system through the Velcro® ends (fig.M1) and then hooking them over the front velcraable band (fig.M2).

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