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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

REF.23770 Ribbed elastic orthosis

REF.23870 Ribbed elastic orthosis with tie-roads

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Item Size	REF.23770B (beige) - REF.23770S (black)					
	XS	S	M	L	XL	XXL
Waist circumference cm	50/60	60/75	75/90	90/105	105/120	120/135
Front height cm	20	20	20	20	20	20
Back height cm	26	26	26	26	26	26

Colour: beige / black

Item Size	REF.23870B (beige) - REF.23870S (black)					
	XS	S	M	L	XL	XXL
Waist circumference cm	50/60	60/75	75/90	90/105	105/120	120/135
Front height cm	20	20	20	20	20	20
Back height cm	26	26	26	26	26	26

Colour: beige / black

MAINTENANCE

⚠ Do not bleach ⚠ No chemical cleaning

⚠ Do not iron ⚠ Do not tumble-dry

🧼 Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

REF.23770B/REF.23770S

- Paravertebral muscle contractures

- Mild traumas

- Mild lumbago

REF.23870B/REF.23870S

- Paravertebral muscle contractures

- Mild traumas

- Mild lumbago

- Mild sciatic lumbalgia and cruralgia

- Post-operative recovery

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Structure made of ribbed elastomer fabric
- 4 lightweight and flexible rods on the back (width 20mm, height 20/22cm)
- Front stays fixed in an oblique position
- Front Velcro® closure

ONLY REF.23770B/23770S

- Reinforced elastic back

ONLY REF.23870B/23870S

- Criss-crossed external elastic tie-rods as a reinforcement, with a Velcro® fastening system

APPLICATION

REF.23770B/ REF.23770S

1 Open the orthosis and place it so that the label is facing upwards

2 Put on the orthosis tensioning it and securing it with the central Velcro® fastener. Make sure the closure is well centered (fig. A).

REF.23870B/ REF.23870S

1 Open the orthosis and temporarily fix the Velcro® straps on the side bands (fig. B).

2 Position the orthosis so that the inner label is facing upwards

3 Put on the orthosis tensioning it and securing it with the central Velcro® fastener. Make sure the closure is well centered (fig. C).

4 Tighten the orthosis by hooking the straps onto the front velcra band.

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Velcro® è un marchio di fabbrica registrato di Velcro Industries B.V.