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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

REF.23001: Postoperative shoe  
REF.23002: Forefoot offloading postoperative shoe  
REF.23003: Heel offloading postoperative shoe  
REF.23004: Toe cap for Stepper

#### DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

#### TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

#### SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

#### WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as proved by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

#### SELECTION/SIZE

| Code      | REF.23001 Stepper surg |         |         |       |         |         |         |
|-----------|------------------------|---------|---------|-------|---------|---------|---------|
| Size      | 3XS                    | XXS     | XS      | S     | M       | L       | XL      |
| Shoe size | 25-32                  | 33-36,5 | 37-38,5 | 39-41 | 41.5-43 | 43.5-45 | 45.5-47 |
| Colour    | black                  |         |         |       |         |         |         |

Sold individually - fits right or left

| Code      | REF.23002 Stepper digit |             |           |         |           |
|-----------|-------------------------|-------------|-----------|---------|-----------|
| Size      | XS                      | S           | M         | L       | XL        |
| Shoe size | 33.5 - 37               | 37.5 - 39.5 | 40 - 41.5 | 42 - 44 | 44.5 - 47 |
| Colour    | black                   |             |           |         |           |

Sold individually - fits right or left

| Code      | REF.23003 Stepper cal |             |           |         |           |
|-----------|-----------------------|-------------|-----------|---------|-----------|
| Size      | XS                    | S           | M         | L       | XL        |
| Shoe size | 33.5 - 37             | 37.5 - 39.5 | 40 - 41.5 | 42 - 44 | 44.5 - 47 |
| Colour    | black                 |             |           |         |           |

Sold individually - fits right or left

| Code   | REF.23004 Stepper toe cap |   |   |   |    |
|--------|---------------------------|---|---|---|----|
| Size   | XS                        | S | M | L | XL |
| Colour | dark grey                 |   |   |   |    |

#### CARE AND MAINTENANCE

- Clean with a sponge dipped in lukewarm water (max. 30°C) and mild soap.
- Dispose of this product and its parts responsibly.

#### INDICATIONS

- Only Stepper surg: For feet requiring stability and offloading after surgery (e.g. osteotomy of the foot, hallux valgus, etc.)
- Only Stepper digit: A postoperative shoe designed to offload pressure from the forefoot.
- Only Stepper cal: Postoperative wear after surgery involving soft tissue and heel bones.
- Only Stepper cal: Targeted offloading relieves pressure following fractures
- Only Stepper cal: Aids the healing of ulcers in the heel area.

#### CONTRAINDICATIONS

There are no known contraindications at this time.

#### CHARACTERISTICS AND MATERIALS

##### SOLE

- The square toe design acts as a bumper to protect the toes.
  - Nonslip sole.
  - The removeable EVA insole can be adapted or replaced if necessary.
  - Only stepper surg: The stiff sole reduces pressure on the heel and forefoot.
  - Only stepper surg: The rocker-style sole facilitates rollover when walking.
  - Only Stepper digit: A 15° wedge offloads pressure from the forefoot.
  - Only Stepper cal: Targeted offloading shifts pressure away from the heel
- ##### UPPER
- Lightweight and breathable.
  - Heel counter for enhanced stability in the heel area.
  - Opens fully so is always quick and easy to don.
  - Pull-back Velcro® strap across the instep with square D-ring and protective padding.
  - Buckle-free symmetrical fastening over the forefoot.
  - Only Stepper surg: Low-profile design so no height compensation is required for the opposite leg.

#### DONNING THE SHOE

- Undo all the Velcro® straps on the shoe.
- Loosen the opening at the front and slide the foot into the shoe.
- Fasten the Velcro® strap across the instep.
- Fasten the wide Velcro® strap over the forefoot.
- The wearer can choose the best way to overlap the tongues in the shoe.

#### HOW TO INSERT THE OPTIONAL TOE CAP REF. 23004

- Remove the insole from inside the postoperative shoe.
- Detach the optional toecap from the Velcro®, remove the protective paper from the sticker and stick it inside the sole, ensuring that it is positioned parallel to the front edge of the sole and no farther than 2 cm from the edge of the sole.
- Position the optional toecap in the front of the shoe and then attach it to the Velcro®.
- Reposition the insole inside the shoe.

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Velcro® è un marchio di fabbrica registrato di Velcro Industries B.V.